



## DINACLICK and Neo-Snap

GE's universal connection systems for both single-tube and dual-tube CRITIKON blood pressure cuffs.



### DINACLICK\*

Innovative and easy-to-use design allows clinicians to quickly connect/disconnect adult and pediatric blood pressure cuffs and provides standardization of your monitors. The audible click lets users know there is a proper connection, and the ergonomically-designed system makes the DINACLICK nearly effortless to use.

#### DINACLICK connector benefits

- **Designed for safety.** May reduce the opportunity for misconnections with luer-style connectors
- **Improved workflow.** Decreases manual connection time
- **Universal connectivity.** Standardization solution for single-tube and dual-tube BP cuffs
- **Intuitive design.** Easy-to-use
- **Non-metallic components.** Material design allows for use in MR environment
- **Patient safety.** Not made with natural rubber latex or DEHP

### Neo-Snap\*

The Neo-Snap connector for neonatal blood pressure cuffs is designed to standards that help prevent unintended disconnect. The Neo-Snap connector provides a secure connection with audible click and tactile feedback.

#### Neo-Snap connector benefits

- **Designed for safety.** Connector meets international standard for improved patient safety
- **Secure connection.** Prevents accidental disconnection
- **Easy-to-use.** Smooth, 360 degree connection that does not have to be aligned to make a secure connection
- **Patient safety.** Not made with natural rubber latex or DEHP

# DINACLICK & Neo-Snap cuff codes

Size	Limb circumference (cm)	Part #	Description
Infant	8 – 13	SFT-P1-2A	SOFT-CUF* INFANT 2T DINACLICK – 20 CUFFS/BOX
Child	12 – 19	SFT-P2-2A	SOFT-CUF CHILD 2T DINACLICK – 20 CUFFS/BOX
Child Long	12 – 19	SFT-P2-2A-L	SOFT-CUF CHILD LONG 2T DINACLICK – 20 CUFFS/BOX
Small Adult	17 – 25	SFT-A1-2A	SOFT-CUF SMALL ADULT 2T DINACLICK – 20 CUFFS/BOX
Small Adult Long	17 – 25	SFT-A1-2A-L	SOFT-CUF SMALL ADULT LONG 2T DINACLICK – 20 CUFFS/BOX
Adult	23 – 33	SFT-A2-2A	SOFT-CUF ADULT 2T DINACLICK – 20 CUFFS/BOX
Adult Long	23 – 33	SFT-A2-2A-L	SOFT-CUF ADULT LONG 2T DINACLICK – 20 CUFFS/BOX
Large Adult	31 – 40	SFT-A3-2A	SOFT-CUF LARGE ADULT 2T DINACLICK – 20 CUFFS/BOX
Large Adult Long	31 – 40	SFT-A3-2A-L	SOFT-CUF LARGE ADULT LONG 2T DINACLICK – 20 CUFFS/BOX
Thigh	38 – 50	SFT-T1-2A	SOFT-CUF THIGH 2T DINACLICK – 20 CUFFS/BOX
5 Small Adult 10 Adult 5 Large Adult		2059303-001	KIT SOFT-CUF 2T DINACLICK – 20 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Adult	26 – 36	SFT-F1-2A	RADIAL-CUF ADULT FOREARM 2T DINACLICK - 20 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Infant	8 – 13	CLA-P1-2A	CLASSIC-CUF* INFANT 2T DINACLICK – 20 CUFFS/BOX
Child	12 – 19	CLA-P2-2A	CLASSIC-CUF CHILD 2T DINACLICK – 20 CUFFS/BOX
Small Adult	17 – 25	CLA-A1-2A	CLASSIC-CUF SMALL ADULT 2T DINACLICK – 20 CUFFS/BOX
Small Adult Long	17 – 25	CLA-A1-2A-L	CLASSIC-CUF SMALL ADULT LONG 2T DINACLICK – 20 CUFFS/BOX
Adult	23 – 33	CLA-A2-2A	CLASSIC-CUF ADULT 2T DINACLICK – 20 CUFFS/BOX
Adult Long	23 – 33	CLA-A2-2A-L	CLASSIC-CUF ADULT LONG 2T DINACLICK – 20 CUFFS/BOX
Large Adult	31 – 40	CLA-A3-2A	CLASSIC-CUF LARGE ADULT 2T DINACLICK – 20 CUFFS/BOX
Large Adult Long	31 – 40	CLA-A3-2A-L	CLASSIC-CUF LARGE ADULT LONG 2T DINACLICK – 20 CUFFS/BOX
Thigh	38 – 50	CLA-T1-2A	CLASSIC-CUF THIGH 2T DINACLICK – 20 CUFFS/BOX
5 Small Adult 10 Adult 5 Large Adult		2059304-001	KIT CLASSIC-CUF 2T DINACLICK – 20 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Neo 1	3 – 6	SFT-N1-1B	SOFT-CUF NEONATE1 1T NEO-SNAP – 20 CUFFS/BOX
Neo 2	4 – 8	SFT-N2-1B	SOFT-CUF NEONATE2 1T NEO-SNAP – 20 CUFFS/BOX
Neo 3	6 – 11	SFT-N3-1B	SOFT-CUF NEONATE3 1T NEO-SNAP – 20 CUFFS/BOX
Neo 4	7 – 13	SFT-N4-1B	SOFT-CUF NEONATE4 1T NEO-SNAP – 20 CUFFS/BOX
Neo 5	8 – 15	SFT-N5-1B	SOFT-CUF NEONATE5 1T NEO-SNAP – 20 CUFFS/BOX
2 – Neo 1 3 – Neo 2 5 – Neo 3 5 – Neo 4 5 – Neo 5		2059305-001	KIT SOFT-CUF NEONATE 1T NEO-SNAP – 20 CUFFS/BOX

# DINACLICK & Neo-Snap cuff codes

Size	Limb circumference (cm)	Part #	Description
Neo 1	3 – 6	SFT-N1-2B	SOFT-CUF NEONATE1 2T NEO-SNAP – 20 CUFFS/BOX
Neo 2	4 – 8	SFT-N2-2B	SOFT-CUF NEONATE2 2T NEO-SNAP – 20 CUFFS/BOX
Neo 3	6 – 11	SFT-N3-2B	SOFT-CUF NEONATE3 2T NEO-SNAP – 20 CUFFS/BOX
Neo 4	7 – 13	SFT-N4-2B	SOFT-CUF NEONATE4 2T NEO-SNAP – 20 CUFFS/BOX
Neo 5	8 – 15	SFT-N5-2B	SOFT-CUF NEONATE5 2T NEO-SNAP – 20 CUFFS/BOX
2 – Neo 1 3 – Neo 2 5 – Neo 3 5 – Neo 4 5 – Neo 5		2059305-002	KIT SOFT-CUF NEONATE 2T NEO-SNAP – 20 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Neo 1	3 – 6	CLA-N1-1B	CLASSIC-CUF NEONATE1 1T NEO-SNAP – 20 CUFFS/BOX
Neo 2	4 – 8	CLA-N2-1B	CLASSIC-CUF NEONATE2 1T NEO-SNAP – 20 CUFFS/BOX
Neo 3	6 – 11	CLA-N3-1B	CLASSIC-CUF NEONATE3 1T NEO-SNAP – 20 CUFFS/BOX
Neo 4	7 – 13	CLA-N4-1B	CLASSIC-CUF NEONATE4 1T NEO-SNAP – 20 CUFFS/BOX
Neo 5	8 – 15	CLA-N5-1B	CLASSIC-CUF NEONATE5 1T NEO-SNAP – 20 CUFFS/BOX
2 – Neo 1 3 – Neo 2 5 – Neo 3 5 – Neo 4 5 – Neo 5		2059306-001	KIT CLASSIC-CUF NEONATE 1T NEO-SNAP – 20 CUFFS/BOX
Neo 1	3 – 6	CLA-N1-2B	CLASSIC-CUF NEONATE1 2T NEO-SNAP – 20 CUFFS/BOX
Neo 2	4 – 8	CLA-N2-2B	CLASSIC-CUF NEONATE2 2T NEO-SNAP – 20 CUFFS/BOX
Neo 3	6 – 11	CLA-N3-2B	CLASSIC-CUF NEONATE3 2T NEO-SNAP – 20 CUFFS/BOX
Neo 4	7 – 13	CLA-N4-2B	CLASSIC-CUF NEONATE4 2T NEO-SNAP – 20 CUFFS/BOX
Neo 5	8 – 15	CLA-N5-2B	CLASSIC-CUF NEONATE5 2T NEO-SNAP – 20 CUFFS/BOX
2 – Neo 1 3 – Neo 2 5 – Neo 3 5 – Neo 4 5 – Neo 5		2059306-002	KIT CLASSIC-CUF NEONATE 2T NEO-SNAP – 20 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Infant	8 – 13	SEN-P1-2A	SENSA-CUF* INFANT 2T DINACLICK – 5 CUFFS/BOX
Child	12 – 19	SEN-P2-2A	SENSA-CUF CHILD 2T DINACLICK – 5 CUFFS/BOX
Small Adult	17 – 25	SEN-A1-2A	SENSA-CUF SMALL ADULT 2T DINACLICK – 5 CUFFS/BOX
Small Adult Long	17 – 25	SEN-A1-2A-L	SENSA-CUF SMALL ADULT LONG 2T DINACLICK – 5 CUFFS/BOX
Adult	23 – 33	SEN-A2-2A	SENSA-CUF ADULT 2T DINACLICK – 5 CUFFS/BOX
Adult Long	23 – 33	SEN-A2-2A-L	SENSA-CUF ADULT LONG 2T DINACLICK – 5 CUFFS/BOX
Large Adult	31 – 40	SEN-A3-2A	SENSA-CUF LARGE ADULT 2T DINACLICK – 5 CUFFS/BOX

# DINACLICK & Neo-Snap cuff codes

Size	Limb circumference (cm)	Part #	Description
Large Adult Long	31 – 40	SEN-A3-2A-L	SENSA-CUF LARGE ADULT LONG 2T DINACLICK – 5 CUFFS/BOX
Thigh	38 – 50	SEN-T1-2A	SENSA-CUF THIGH 2T DINACLICK – 5 CUFFS/BOX
1 Small Adult 1 Adult 1 Large Adult		2059300-001	KIT SENSA-CUF 2T DINACLICK – 3 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Infant	8 – 13	DUR-P1-2A	DURA-CUF* INFANT 2T DINACLICK – 5 CUFFS/BOX
Child	12 – 19	DUR-P2-2A	DURA-CUF CHILD 2T DINACLICK – 5 CUFFS/BOX
Child Long	12 – 19	DUR-P2-2A-L	DURA-CUF CHILD LONG 2T DINACLICK – 5 CUFFS/BOX
Small Adult	17 – 25	DUR-A1-2A	DURA-CUF SMALL ADULT 2T DINACLICK – 5 CUFFS/BOX
Small Adult Long	17 – 25	DUR-A1-2A-L	DURA-CUF SMALL ADULT LONG 2T DINACLICK – 5 CUFFS/BOX
Adult	23 – 33	DUR-A2-2A	DURA-CUF ADULT 2T DINACLICK – 5 CUFFS/BOX
Adult Long	23 – 33	DUR-A2-2A-L	DURA-CUF ADULT LONG 2T DINACLICK – 5 CUFFS/BOX
Large Adult	31 – 40	DUR-A3-2A	DURA-CUF LARGE ADULT 2T DINACLICK – 5 CUFFS/BOX
Large Adult Long	31 – 40	DUR-A3-2A-L	DURA-CUF LARGE ADULT LONG 2T DINACLICK – 5 CUFFS/BOX
Thigh	38 – 50	DUR-T1-2A	DURA-CUF THIGH 2T DINACLICK – 5 CUFFS/BOX
Infant	8 – 13	DUR-P1-2A	DURA-CUF INFANT 2T DINACLICK – 1 CUFF/PACK
Child	12 – 19	DUR-P2-2A	DURA-CUF CHILD 2T DINACLICK – 1 CUFF/PACK
Child Long	12 – 19	DUR-P2-2A-L	DURA-CUF CHILD LONG 2T DINACLICK – 1 CUFF/PACK
Small Adult	17 – 25	DUR-A1-2A	DURA-CUF SMALL ADULT 2T DINACLICK – 1 CUFF/PACK
Small Adult Long	17 – 25	DUR-A1-2A-L	DURA-CUF SMALL ADULT LONG 2T DINACLICK – 1 CUFF/PACK
Adult	23 – 33	DUR-A2-2A	DURA-CUF ADULT 2T DINACLICK – 1 CUFF/PACK
Adult Long	23 – 33	DUR-A2-2A-L	DURA-CUF ADULT LONG 2T DINACLICK – 1 CUFF/PACK
Large Adult	31 – 40	DUR-A3-2A	DURA-CUF LARGE ADULT 2T DINACLICK – 1 CUFF/PACK
Large Adult Long	31 – 40	DUR-A3-2A-L	DURA-CUF LARGE ADULT LONG 2T DINACLICK – 1 CUFF/PACK
Thigh	38 – 50	DUR-T1-2A	DURA-CUF THIGH 2T DINACLICK – 1 CUFF/PACK
1 Small Adult 1 Adult 1 Large Adult		2059301-001	KIT DURA-CUF 2T DINACLICK – 3 CUFFS/BOX

Size	Limb circumference (cm)	Part #	Description
Infant	8 – 13	CLA-P1-2A-X	CLASSIC-CUF INFANT 2T DINACLICK ISO – 20 CUFFS/BOX
Child	12 – 19	CLA-P2-2A-X	CLASSIC-CUF CHILD 2T DINACLICK ISO – 20 CUFFS/BOX
Small Adult	17 – 25	CLA-A1-2A-X	CLASSIC-CUF SMALL ADULT 2T DINACLICK ISO – 20 CUFFS/BOX
Adult	23 – 33	CLA-A2-2A-X	CLASSIC-CUF ADULT 2T DINACLICK ISO – 20 CUFFS/BOX
Large Adult	31 – 40	CLA-A3-2A-X	CLASSIC-CUF LARGE ADULT 2T DINACLICK ISO – 20 CUFFS/BOX
Thigh	38 – 50	CLA-T1-2A-X	CLASSIC-CUF THIGH 2T DINACLICK ISO – 20 CUFFS/BOX

# Adapters and accessories

Part #	Adapters for air hoses	Quantity/package
2059291-001	Y-ADAPTER: 1T MALE BAYONET TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059291-002	Y-ADAPTER: 1T FEMALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059291-003	Y-ADAPTER: 1T MALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059291-004	Y-ADAPTER: 1T UNFITTED TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059292-001	X-ADAPTER: 2T MALE BAYONET TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059292-002	X-ADAPTER: 2T FEMALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059292-003	X-ADAPTER: 2T MATED SUBMIN TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059292-004	X-ADAPTER: 2T MALE SCREW TO 2T DINACLICK FEMALE RECEPT	10 ADAPTERS/PACK
2059290-001	Y-ADAPTER: 1T MALE SLIP LUER TO 2T NEO-SNAP MALE PLUG	10 ADAPTERS/PACK
2059290-002	Y-ADAPTER: 1T MALE MICRO PLUG TO 2T NEO-SNAP MALE PLUG	10 ADAPTERS/PACK
2059297-001	ADAPTER: 1T MALE MICRO PLUG TO 1T NEO-SNAP MALE PLUG	10 ADAPTERS/PACK

Part #	Wall manometer and adapters for manual inflation systems	Quantity/package
2059289-001	ADAPTER: WALL MANOMETER W/ COIL TUBING, BULB AND VALVE TO 2T DINACLICK FEMALE RECEPT	5 ADAPTERS/BOX
2059289-002	ADAPTER: WALL MANOMETER W/ 13" TUBING TO 2T DINACLICK FEMALE RECEPT	5 ADAPTERS/PACK
2059298-001	WALL MANOMETER W/ CUFF	1 WALL MANOMETER & 1 CUFF/BOX
2059299-001	WALL MANOMETER W/ CUFF & ROLLSTAND	1 WALL MANOMETER, 1 CUFF & 1 ROLLSTAND/BOX

Part #	Accessories	Quantity/package
2059293-001	2T NEO-SNAP MALE PLUG	10 PLUGS/BOX
2059294-001	2T DINACLICK FEMALE RECEPT	10 DINACLICK/PACK
2059295-001	WASH PLUGS - FOR NEO-SNAP FEMALE & DINACLICK MALE	50 PLUGS/BOX
2055630-001	1T DINACLICK 1/8" TUBE FEMALE RECEPT	10 DINACLICK/PACK
2055630-002	1T CONN 2-1 DINACLICK 5/32" TUBE FEMALE RECEPT	10 DINACLICK/PACK

# Hoses

Hose type	Length	Hose	Description	Connector style
Adult	1.2 M/4 FT	2058203-001	CARE HOSE ADULT 2T DINACLICK	2T DINACLICK
Adult	3.6 M/12 FT	2058203-002	CARE HOSE ADULT 2T DINACLICK	2T DINACLICK
Adult	7.2 M/24 FT	2058203-003	CARE HOSE ADULT 2T DINACLICK	2T DINACLICK
Adult	3.6 M/12 FT	2058205-001	CRITIKON* HOSE ADULT 2T DINACLICK	2T DINACLICK
Adult	7.2 M/24 FT	2058205-002	CRITIKON HOSE ADULT 2T DINACLICK	2T DINACLICK
Neo	1.2 M/4 FT	2058204-001	CARE HOSE NEO 2T NEO-SNAP	2T NEO-SNAP
Neo	2.4 M/8 FT	2058204-002	CARE HOSE NEO 2T NEO-SNAP	2T NEO-SNAP
Neo	3.6 M/12 FT	2058204-003	CARE HOSE NEO 2T NEO-SNAP	2T NEO-SNAP
Neo	6.0 M/19.7 FT	2058204-004	CARE HOSE NEO 2T NEO-SNAP	2T NEO-SNAP
Neo	3.6 M/12 FT	2058206-001	CRITIKON HOSE NEO 2T NEO-SNAP	2T NEO-SNAP

## Imagination at work

GE Healthcare  
8200 West Tower Avenue  
Milwaukee, WI 53223  
U.S.A.

[cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information. © 2014 General Electric Company - All rights reserved. GE and GE Monogram are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. \*Trademark of General Electric Company.

DOC1462023 NON-EU 02/14