



CRITIKON CLASSIC-CUF

Blood Pressure Cuffs

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong to withstand repeated inflations – yet economical enough for single patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large, bold print helps you quickly choose the right cuff, enabling more accurate measurement

Constructed for patient safety and durability

- Easy-to-clean, vinyl material does not contain natural rubber latex or DEHP, making them safer for your patients
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility

Cuff assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



CLASSIC-CUF			ISO 80369-5 compliant	
				
Size	Color	Limb Circumference	2-Tube DINACLICK Connector	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	CLA-P1-2A-5S	2300
Child	Green	12-19 cm	CLA-P2-2A-5S	2301
Small Adult	Light Blue	17-25 cm	CLA-A1-2A-5S	2302
Adult	Navy	23-33 cm	CLA-A2-2A-5S	2303
Adult Long	Navy	23-33 cm	CLA-A2-2A-L-5S	2306
Large Adult	Rose	31-40 cm	CLA-A3-2A-5S	2304
Large Adult Long	Rose	31-40 cm	CLA-A3-2A-L-5S	2307
Thigh	Brown	38-50 cm	CLA-T1-2A-5S	2305

Note: 20 cuffs/box

DINACLICK Assortment Packs

Connector Type	Part Number (ISO 80369-5 compliant)	Description
2-Tube DINACLICK Connectors	2059304-501	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2059304-502	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/pack)



© 2020 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CLASSIC-CUF, CRITIKON and DINACLICK are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00921XE 8/20