



Are your labor and delivery transducers reliable when you need it most?

Corometrics™ Nautilus Nine-Crystal transducers help provide reassurance and support

Available directly from GE Healthcare, Corometrics transducers are available for replacement throughout the life of your Corometrics monitor. The combo transducers, for ultrasound and toco, offer a novel combination of both button and loop connection. The transducers are designed to support minimal repositioning even with patients of varying BMI and gestational age.

Not all transducers or transducer repairs are compatible with your Corometrics monitor. Corometrics Nautilus transducers are the only verified transducers for use with Corometrics 250cx and 170 series monitors.

Ultrasound transducer button and loop

Part number: 2108346-001

Toco transducer button and loop

Part number: 2108347-001

Abdominal belt

Part number: 2015827-001; 10/Pkg

Part number: 2015919-001; 200/Pkg

To order:

- Visiting www.services.gehealthcare.co.uk
- Email gehccaccessorysales@ge.com with PO number
- Calling Customer Service at 0800-032-9201, option 2, then 3

gehealthcare.co.uk/products/accessories-and-supplies



Patient Comfort

The flat toco transducer is designed for increased comfort and is watertight for cleanability and in-water use.*



Parameter Excellence

Nine-crystal technology helps achieve a wide focal region for more uniform coverage at greater depths.

© 2021 General Electric Company – All rights reserved.

*When used with Mini Telemetry

Product color may vary from image shown.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Corometrics are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company.

JB01092XE Rev 2 01/21