



# CRITIKON DURA-CUF

## Blood Pressure Cuffs

CRITIKON™ DURA-CUF™ is gentle and designed for durability. Made with tough, flexible nylon, they are convenient and easy-to-clean with curved corners to enhance your patient's comfort. They are constructed for multiple patient use in low-acuity care settings and ideal for physicians offices where the risk of contamination is diminished.



### Designed for patient safety, comfort and accuracy

- Made with durable materials that can withstand extensive use, and multiple disinfections and cleanings
- Designed with rounded corners to enhance patient comfort
- Constructed with materials not made with PVC or natural rubber latex, making them safer for your patient
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment

### Suitable for virtually every physician's office

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors
- A full range of adapters and connectors are compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility
- Three types of cuff assortment packs enable easier ordering and optimization of inventory space



# DURA-CUF

ISO 80369-5 compliant



Size	Color	Limb Circumference	2-Tube DINACLICK Connector	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	DUR-P1-2A-5S	002274
Child	Green	12-19 cm	DUR-P2-2A-5S	002275
Child Long	Green	12-19 cm	DUR-P2-2A-L-5S	
Small Adult	Royal Blue	17-25 cm	DUR-A1-2A-5S	002276
Small Adult Long	Royal Blue	17-25 cm	DUR A1-2A-L-5S	
Adult	Navy	23-33 cm	DUR-A2-2A-5S	002277
Adult Long	Navy	23-33 cm	DUR-A2-2A-L-5S	002130
Large Adult	Wine	31-40 cm	DUR-A3-2A-5S	002278
Large Adult Long	Wine	31-40 cm	DUR-A3-2A-L-5S	002131
Thigh	Brown	38-50 cm	DUR-T1-2A-5S	002279

Note: 5 cuffs/box

## DINACLICK Assortment Packs

Connector Type	Part Number (ISO 80369-5 compliant)	Description
2-Tube DINACLICK Connectors	2059301-501	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)
	2059301-502	2 Infant, 3 Child, 1 Small Adult (6 cuffs/box)
	2059301-503	1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult (6 cuffs/box)
	2059301-504	1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh (6 cuffs/box)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [gehealthcare.co.uk/products/patient-monitoring/patient-monitors/non-invasive-blood-pressure-cuffs](https://gehealthcare.co.uk/products/patient-monitoring/patient-monitors/non-invasive-blood-pressure-cuffs)



© 2020 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, CRITIKON, DURA-CUF and DINACLICK are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB00923XE 8/20