


Versana Premier™

Easy to use. Easy to learn. Easy to diagnose.



Some innovations feel so natural from the start,
you wonder how you ever lived without them



Versana Premier

Easy to use. Easy to learn. Easy to diagnose.


Investing in and introducing a new ultrasound system to your practice can be a significant and challenging decision. Versana Premier helps clinicians enjoy the sophisticated yet familiar simplicity of an ultrasound system that feels natural from the start, making it easier than ever to deliver quality care to a wide range of patients every day.


Easy to use and easy to learn, Versana Premier lets clinicians easily and confidently diagnose a wide spectrum of orthopedic conditions. Whether you're in an orthopedic clinic, independent imaging center, hospital radiology department, general practice office, health checkup center or other facility, you can count on Versana Premier as a reliable and user-friendly ultrasound system designed to empower clinicians, enhance their efficiency, and provide optimal care, all with effortless productivity in mind.


Easy to use

From its design to its automated tools to its reliable service and support, Versana Premier delivers an effortless user experience.


Transform your workflow and accelerate scanning with advanced, intuitive features

 Experience up to 38% faster workflow with *Whizz clinical features*, saving time and simplifying workflow¹

 Create your own report templates with *Whizz Report* and import and display clinical guidelines for references while scanning with *Whizz Notes*

 Use new on-console features, such as *Assistant*, to fit seamlessly into your daily routine

Advanced design and capabilities to simplify operations and meet the needs of your practice and patients

 Easily move the system around your facility with its portable design that is equipped with enhanced ergonomic updates, improved accessory management, and extended battery-powered scanning



Easy to learn

Versana Premier offers an array of clinical and product education resources to help you and your staff master the system quickly, broaden your range of exam capabilities, and spend more quality time with your patients.

Robust educational support, advanced tools, & reliable service provide a comprehensive solution for users' practices both now & in the future



Completely backed by GE HealthCare for education and service with *Versana Club*, *My Trainer*, and *eDelivery*



1 — Maximize your efficiency and train new users with built-in reference and guidance tools, including *Scan Assistant*, *Scan Coach* now for *MSK*, and *Scan Coach Voice Comments*



My Trainer on-board help walks you through the setup, operation, cleaning and maintenance of your Versana ultrasound system. It helps you get started with your Versana system quickly



Receive real-time reference information to help locate the desired scan plane with *Scan Coach now for MSK*, including 3D animations and illustrations showing the target anatomy and suggested probe placement and reference images showing how the scanned image may appear on screen



Easy to diagnose

Powered by *VisionBoost Architecture*, Versana Premier is designed to provide optimal imaging and tools needed to perform broad range of MSK ultrasound exams, so that users can accurately make decisions and reach a diagnosis with the greatest of ease.

Capture exceptional images with advanced clinical features to address challenging MSK cases



Optimize images while scanning with *Whizz* clinical features and easily improve Color Doppler with *Whizz Color Flow*



Get additional clinical information about the relative tissue stiffness with *Strain Elastography*



Shift between Smooth Mode and Sharp Mode to view image features as preferred with *Whizz Easy Style*



Assess for inflammation or soft tissue neoplasms with *B-Flow, B-Flow Color* and *Power Doppler Imaging*



Leverage a wide range of transducers, including new *L4-20t-RS*, which produces high resolution images across clinical specialties and care areas

Enhance your visualization and track image changes to help improve diagnostic confidence



Virtual Convex enables a larger field of view in the far field while *LOGIQView* gives the possibility to do panoramic imaging by moving the probe longitudinal along a structure



Visualize, evaluate and take measurements on the full screen for selected ROI with *V-Zoom*



Adjust needle gain and angle when you add *Needle Recognition Software*, which provides precise and accurate display of the needle, anatomy and motion




Easily compare current and prior exams and track treatment progress with *Follow-Up Tool*



Effortlessly adapt to various orthopedic conditions with a wide range of probes*

*Probes may not be available in all countries and regions

Ergonomics & hardware updates




*Decibel range of 21.96 to 39.51.
Equivalent to Versana Premier V2

Versana Premier


is an intelligently productive, user-friendly ultrasound system with multiple exam protocols, tools, and advanced applications to streamline processes so it's easier than ever for you to deliver quality care to your wide range of patients every day.




Image gallery




L4-20t-RS Median Nerve




L4-20t-RS Patellar Tendon with Virtual Convex




L8-18i-RS Median Nerve with Follow-Up Fusion



12L-RS Supraspinatus Tendon



12L-RS Biceps Tendon



L3-12-RS Wrist 1st Compartment


Familiar from day one

An ultrasound solution making the delivery of quality patient care easier than ever before

Advance your ultrasound skills with Versana Club

No matter your skill level, you can learn to harness the full power of your Versana Premier system with a membership in the Versana Club. The Versana Club offers resources to help you build your skills, network with others and share your ultrasound knowledge. Membership perks include:

- Image-sharing with peers around the world
- Web-based video tutorials and other educational offerings
- Information on Versana service solutions and products
- Articles on ultrasound and patient care



Register now at
www.versanaclub.net/emea/home

Customized solutions and path to growth

Versana Premier can scale to suit your practice today and grow with you. Select a basic console and software. Then choose from a variety of optional upgrades as your practice advances

Experience the strength of the GE HealthCare family

Your support team is fully trained, GE HealthCare qualified, deeply experienced, and nearby.

- Expert support during and after purchase
- In-depth staff product training



Versana Essential™

Versana Premier™

Versana Balance™

Versana Active™





1. GE HealthCare-sponsored usability study, using Whizz. Study was conducted in December 2022. JB24263XX

© 2024 GE HealthCare. Versana Active, Versana Essential, Versana Balance, Versana Premier are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license. JB01361AFG
Commercial availability of GE HealthCare medical systems is subject to meeting local requirements in a given country or region. Contact a GE HealthCare Representative for more information.
Intended for healthcare professionals only.
August 2024



GE HealthCare