

CRITIKON CLASSIC-CUF

Blood Pressure Cuffs

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong to withstand repeated inflations – yet economical enough for single patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE HealthCare patient monitors, with repeated inflations to ensure accurate readings.
- Large, bold print helps you quickly choose the right cuff, enabling more accurate measurement.

Constructed for patient safety and durability

- Easy-to-clean, vinyl material does not contain natural rubber latex or DEHP, making them safer for your patients.
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit.
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment.

Designed to improve workflow

- DINACLICK's easy-to-use design allows clinicians to quickly connect/ disconnect adult and pediatric blood pressure cuffs.
- Compliant with AAMI standard and AHA sizing guidelines.
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility.

Cuff assortment packs available

- Assortment packs are available to be flexible for all users and simplified ordering.



CLASSIC-CUF			ISO 80369-5 compliant	
				
Size	Color	Limb Circumference	2-Tube DINACLICK Connector	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	CLA-P1-2A-5S	2300
Child	Green	12-19 cm	CLA-P2-2A-5S	2301
Small Adult	Light Blue	17-25 cm	CLA-A1-2A-5S	2302
Adult	Navy	23-33 cm	CLA-A2-2A-5S	2303
Adult Long	Navy	23-33 cm	CLA-A2-2A-L-5S	2306
Large Adult	Rose	31-40 cm	CLA-A3-2A-5S	2304
Large Adult Long	Rose	31-40 cm	CLA-A3-2A-L-5S	2307
Thigh	Brown	38-50 cm	CLA-T1-2A-5S	2305

Note: 20 cuffs/box

DINACLICK Assortment Packs

Connector Type	Part Number (ISO 80369-5 compliant)	Description
2-Tube DINACLICK Connectors	2059304-501	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2059304-502	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/pack)