

# CRITIKON Neonatal Blood Pressure Cuffs

Providing big benefits for your tiniest patients, GE HealthCare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl. The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single patient use.



## A trusted choice for clinical excellence

- Tested for performance accuracy and clinically validated to help you deliver quality patient care
- A full range of neonatal sizes in playful, color-coded animal designs to help you find the right size at a glance, enabling more accurate measurement
- Choice of extra soft, absorbent material (SOFT-CUF™) or easy-to-clean vinyl (CLASSIC-CUF™)

## Constructed for patient safety

- Single-patient-use for pairing with a single baby throughout their hospital stay to help with your infection control efforts
- Cuff materials are not made with natural rubber latex or DEHP, making them safer for the babies
- Neo-Snap™ connector meets international standards, prevents misconnections to luer fittings, and helps prevent accidental disconnections for patient safety



Neo-Snap Connector

## Designed to improve workflow

- Neo-Snap connector does not have to be aligned to make a secure connection, enabling you to focus more on the baby
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors

## Cuff assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



Neonatal			2-Tube Neo-Snap Connector		1-Tube Neo-Snap Connector	
			SOFT-CUF	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF
Size	Color	Limb Circumference				
1	Orange	3-6 cm	SFT-N1-2B	CLA-N1-2B	SFT-N1-1B	CLA-N1-1B
2	Light Blue	4-8 cm	SFT-N2-2B	CLA-N2-2B	SFT-N2-1B	CLA-N2-1B
3	Green	6-11 cm	SFT-N3-2B	CLA-N3-2B	SFT-N3-1B	CLA-N3-1B
4	Navy	7-13 cm	SFT-N4-2B	CLA-N4-2B	SFT-N4-1B	CLA-N4-1B
5	Rose	8-15 cm	SFT-N5-2B	CLA-N5-2B	SFT-N5-1B	CLA-N5-1B

Note: 20 cuffs/box

Neonatal Assortment Pack	Cuff Type	Part Number	Description
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-003	1 Neo #3, 1 Neo #4, 1 Neo #5 (3 cuffs/pack)
1-Tube Neo-Snap Connectors	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)

To learn more about GE HealthCare comprehensive cuff solutions, visit us at [gehealthcare.co.uk/products/clinical-accessories/non-invasive-blood-pressure-cuffs](https://gehealthcare.co.uk/products/clinical-accessories/non-invasive-blood-pressure-cuffs)

© 2024 GE HealthCare

GE HealthCare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE HealthCare representative for the most current information. CRITIKON, SOFT-CUT and DINACLICK are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

JB23415XX

