



GE Healthcare

Reusable Breathing Bags



Reusable breathing bags from GE Healthcare are tested and verified to function with GE Healthcare anesthesia machines when applicable. The products perform optimally and are produced from high quality silicone materials to ensure product durability and usability.

FEATURES

- Reusable silicone bags are autoclavable up to 50 times with high quality materials and feel
- All breathing bags use standard 22mm cuff for ease of assembly
- Available in a range of sizes to suit both patient size and clinical user preferences
- Bag loop included on all sizes for ease of use and storage

**Our equipment.
Our accessories.
Simply compatible.**

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

REUSABLE BREATHING BAG OFFERINGS

Silicone reusable breathing circuit bags are used in operating rooms and connect to GE Healthcare anesthesia machines to support the patient with their respiration when the operator has selected modes for manual ventilation.

PRODUCT SPECIFICATIONS

Material composition

Breathing Bags Silicone

Latex content

Breathing Circuit Bags do not contain natural rubber latex

DEHP content

Does not contain phthalate DEHP

Product use

Breathing Circuit Bags are reusable devices. They can be reused up to 50 times.

Sterility

Breathing Circuit Bags are supplied non-sterile

Storage

Store in a cool, dry place out of direct sunlight

Shelf life

Shelf life of 5 years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

PART LIST

Part Number	Description	Quantity
2106570-001	Manual bag, 0.5L, reusable, 22mm cuff	1
2106570-002	Manual bag, 1.0L, reusable, 22mm cuff	1
2106570-003	Manual bag, 2.0L, reusable, 22mm cuff	1
2106570-004	Manual bag, 3.0L, reusable, 22mm cuff	1
2106570-025	Manual bag, 0.5L, reusable, 22mm cuff	5
2106570-026	Manual bag, 1.0L, reusable, 22mm cuff	5
2106570-027	Manual bag, 2.0L, reusable, 22mm cuff	5
2106570-028	Manual bag, 3.0L, reusable, 22mm cuff	5

Disposal considerations


Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Packaging materials

Primary	Polybag – Polyethylene (PE)
Secondary	Carton / Box – Cardboard

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com

Distributed by:
Datex-Ohmeda, Inc.

 *Manufactured by:*
Flexicare Medical Limited
Cynon Valley Business Park,
Mountain Ash, CF45 4ER, UK