



Quick connections

Cardiac Output Trunk Cables

Cardiac output cables from GE Healthcare are tested and verified to function with certain GE Healthcare monitoring equipment. Parameter color-coding makes them quick to connect.

- Cable curls and tangles are minimized by premium quality materials
- A molded connector assembly helps keep cables flexible and prevents ingress of fluids
- The cables have a connector for an injectate probe and also a standard three-pin catheter connector that fits most Swan-Ganz™ catheters

Features

GE cardiac output cables are available in lengths of 3.6 m/12 ft. and 2.4 m/8 ft.

Product compatibility

Cables are compatible and verified with GE CARESCAPE™ B450, B650 and B850 Monitors via the CARESCAPE Patient Data Module (PDM).

Accessories simply compatible with your equipment

You can trust GE Healthcare accessories and consumables to keep your clinical devices working at peak efficiency.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step

PRODUCT SPECIFICATIONS

Materials

Not made with natural rubber latex or DEHP
ROHS (Restriction of hazardous substances) compliant

Operating conditions

Check the host device manual for recommendations.

Storage conditions

Recommended temperature -10° to 50°C (14° to 122°F)

Cleaning solutions

Ethyl alcohol maximum 90% by volume

Isopropyl alcohol maximum 60% by weight

Glutaraldehyde 2%

Phenol 2%

Do not sterilize or autoclave.

Compliance

CE Mark

ORDERING INFORMATION

| Part Number | Description |
|-------------|------------------------------------|
| 2104180-001 | Cardiac Output Cable, 3.6 m/12 ft. |
| 2104180-002 | Cardiac Output Cable, 2.4 m/8 ft. |

Product may not be available in all countries and regions.
Contact a GE Healthcare Representative for more information.
Please visit www.gehealthcare.com/promotional-locations

Imagination at work

© 2018 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram, are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. Swan-Ganz is a trademark of EDWARDS LIFESCIENCES CORPORATION. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2146445 6/18