

Clinical Accessories

Trusted. Reliable. Simple.



Contents

PATIENT MONITORING	4				
ECG	6				
Monitoring Recording Paper	6				
Trunk Cables	7				
Leadwires	8-10				
Trunk Cable Accessories	11				
Pulse Oximetry	12				
Sensors & Cables	12-18				
Temperature	19				
Cables & Probes	19-20				
Invasive Blood Pressure	21				
Cables	21-22				
Cuffs	23				
Hoses	23				
DINACLICK	24				
SOFT-CUF	25				
CLASSIC-CUF	26				
DURA-CUF	27				
YELLOW CLASSIC-CUF	28				
RADIAL-CUF	29				
Neonatal	30				
Cardiac Output	31				
Cables	31-32				
Other Accessories	32				
Adequacy of Anesthesia	33				
Cable	33-34				
Sensors	33-34				
Neurology	35				
EEG Leadsets	35				
EEG Electrodes	36				
AEP Accessories	36				
Patient Spirometry	37				
Sensors	37				
Tubes	38				
Kits	38				
Gas Monitoring	39				
Airway Adapter Sets	39				
Gas Sampling Lines	40				
Airway Adapters	40				
Water Traps	40				
ANESTHESIA & RESPIRATORY CARE	41				
Absorbents	43				
Absorbers	43				
Canisters	43				
Dust Filter	43				
Vaporizer Bottle Adapters	43				
Filters & Breathing Masks	44				
Reservoir Bags & Gas Scavenging	45				
Anesthesia Breathing Bags	45				
Scavenging Bags, Tubes, & Connectors	45				
Patient Breathing Circuits	46				
Circuit Startup Kits	46				
Circuit Kits	47				
Circuits	48-49				
Circuit Components	50				
Flow Sensors – Respiratory	51				
Respiratory Flow Sensors and Expiratory Valve Assemblies	51				
Flow Sensors – Anesthesia	52				
Nebulizer Accessories	53				
Anesthesia O2 Cells	53				
Respiratory Safety Guard	54				
Pressure Catheters	54				
DIAGNOSTIC CARDIOLOGY	55				
MAC™, CASE™ and CardioSoft™	57				
ECG	57				
Cables	57				
Leadwires	57-58				
Accessories	59				
Recording Paper	60				
Electrodes	61				
Holter Monitors	62				
ECG – CM 4000	62				
Leadwire Sets	62				
Accessories	63				
ECG – SEER® 12	63				
Leadwire Sets	63				
Accessories	63				
ECG – SEER® 1000	64				
Leadwire Sets	64				
Accessories	65				
ECG – SEER® Light	66				
Leadwire and Cable Kits	66				
Accessories	67				
MATERNAL INFANT CARE	68				
Maternal Fetal Care	70				
Mini Telemetry	70				
Recording Thermal Paper	70				
Remote Event Marker	70				
Novii™	71				
Novii™ Patch	71				
Cables	71				
Corometrics™	72				
Cables	72				
Transducer Belts	72				
Ultrasound Transducers	72				
Internal Monitoring	73				
Cables	73				
Electrode	73				
Pad	73				
Infant Care	74				
Giraffe™/Panda™	74-75				
Air filters	74				
Clip	74				
Cover	75				
Phototherapy Accessories	76				
Pad Covers	76				
Resuscitation	77				
T-Piece	77				
Lullaby™ Warmer	78				
Temperature Probes	78				
Temperature Probe Patches	78				
Shuttle™	78				
Patient Restraint	78				

Welcome

GE Healthcare offers you a simple, reliable one-stop solution for high-quality GE Healthcare-compatible clinical accessories. We provide an expanded inventory of products that have been verified by our engineers. The compatible devices listed may not be available for sale and/or cleared/approved in all countries. Check with your representative for product availability.

PATIENT MONITORING



PATIENT MONITORING

ECG

- Leadwires
- Recording Paper
- Trunk Cable
- Trunk Cables Accessories

Pulse Oximetry

- Cables
- Clip
- Other Accessories
- Probes
- Sensors
- Tape & Wrap

Temperature

- Cables
- Probes

Invasive Blood Pressure

- Cables

Cuffs

- Hoses
- DINACLICK
- SOFT-CUFF
- DURA-CUF
- CLASSIC-CUF
- YELLOW CLASSIC-CUF
- RADIAL-CUF
- Neonatal

Cardiac Output

- Cables
- Catheter Kits
- Other Accessories

Adequacy of Anesthesia

- Cable
- Electrodes
- Sensors

Neurology

- AEP Accessories
- EEG Electrodes
- EEG Leadsets

Patient Spirometry

- Kits
- Sensors
- Tubes

Gas Monitoring

- Airway Adapter Sets
- Airway Adapters
- Gas Sampling Lines
- Water Traps

ECG



MONITORING RECORDING PAPER



MS Recording Thermal Paper

Length	CARESCAPE™ B450/B650 Monitors, CARESCAPE™ B850 Monitor with PRN 50-M+ printer	B20(i)/B40(i) Monitors	B1x5 Monitors	CARESCAPE™ V100 and VC150 Monitors	Part Number	Quantity
57.2 mm x 24.4 m				X	2106825-001	10
50 mm x 25.6 m	X	X	X		2106824-001	20
50 mm x 27.4 m	X	X	X		2106823-001	48

For ECG paper for Diagnostic cardiology and Maternal infant care equipment, please go to pages 58 and 68.

Not available in all countries



2106305-003



2106310-002

ECG Trunk Cables – IEC			
Description	Length	Part Number	Quantity
Neonatal DIN 3-lead	3.6 m	2106306-003	1
	1.2 m	2106306-004	1
3/5-lead	3.6 m	2106305-003	1
	1.2 m	2106305-004	1
	6 m	2106311-002	1
3/5-lead w/ESU filter	3.6 m	2106310-002	1
6-lead	3.6 m	2106307-003	1
	1.2 m	2106307-004	1
12-lead	3.6 m	2106308-003	1
	1.2 m	2106308-004	1



2106307-003



2106308-004

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module E-PSM
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ Monitor with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry

Specifications available on request



ECG Trunk Cables – IEC			
Description	Length	Part Number	Quantity
Neonatal trunk cable	3.6 m	2080557-003	1
	1.2 m	2080557-004	1

Not available in all countries



ECG Leadwires Grabber – IEC Compliant

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm*	2106390-003	1
		130 cm	2106390-004	
5-lead	Molded	74 cm*	2106389-003	1
		130 cm	2106389-004	
		mix 74 cm, 130 cm	2106389-006	
	Non-molded	74 cm *	2106391-003	
		130 cm	2106391-004	
		mix 74 cm, 130 cm	2106391-006	
5-lead C2-6	Non-molded	130 cm	2106393-004	1
		74 cm	2106393-002	
6-lead	Non-molded	74 cm*	2106398-001	1
		130 cm	2106398-002	
		mix 74 cm, 130 cm	2106398-003	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

*74cm 3-, 5-, and 6-leadwire sets are compatible with ApexPro™

Not available in all countries





2106385-004



2106383-006



2106381-004



2106394-002



2106398-001

ECG Leadwires Snap – IEC Compliant

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm*	2106385-003	1
		130 cm	2106385-004	
5-lead	Molded	74 cm*	2106383-003	1
		130 cm	2106383-004	
		mix 74 c m, 130 cm	2106383-006	
	Non-molded	74 cm*	2106381-003	
130 cm		2106381-004		
5-lead C2-6	Non-molded	130 cm	2106394-004	1
		74 cm	2106394-002	
6-lead	Non-molded	74 cm*	2106396-001	1
		130 cm	2106396-002	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

*74cm 3-, 5-, and 6-leadwire sets are compatible with ApexPro™

Not available in all countries





ECG Leadwires – IEC Compliant				
Description	Connection	Length	Part Number	Quantity
4-lead	Snap	74 cm/29 in	2066468-002	1
		130 cm/51 in	2066468-006	
	Grabber	74 cm/29 in	2066468-004	1
		130 cm/51 in	2066468-008	

Product compatibility

- CARESCAPE™ ONE Monitor

Specifications available on request

Not available in all countries



2106403-003



2106403-001



2106399-001



2106399-002

ECG Trunk Cable Accessories

Description		Part Number	Quantity
DIN Converter	3-Lead	2106403-003	1
	5-Lead	2106403-001	1
Leadwire separator		2106399-001	3
Bedsheet clip		2106399-002	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry

Specifications available on request

Not available in all countries



Pulse Oximetry



TRUSIGNAL™

GE Healthcare provides compatibility with three major brands of Pulse Oximetry (SpO₂) measurement to meet all of your needs:

- GE Healthcare TruSignal™
- MASIMO
- Nellcor



TS-F2-GE



TS-SA4-GE



TS-E2-GE

TruSignal™ Integrated Reusable Sensors						
Description	Connector	Patient	Patient Weight	Length	Part Number	Quantity
Finger	GE	Adult, Pediatric	≥ 20 kg	2 m	TS-F2-GE	1
				4 m	TS-F4-GE	1
	Ohmeda			1 m	TS-F1-H	1
				Datex	4 m	TS-F4-N
FingerTip	GE	Adult	> 30 kg		4 m	TS-SA4-GE
PediTip		Pediatric	15-30 kg	3 m	TS-SP3-GE	1
Ear	GE	Adult, Pediatric	≥ 10 kg	2 m	TS-E2-GE	1
				Ohmeda	4 m	TS-E4-GE
	TS-E4-H					1
	TS-E4-N					1
	Datex					

Not available in all countries





TS-F-D



TS-SA-D/TS-SP-D




TS-E-D



TS-W-D



TS-SE-3



TruSignal™ Reusable Sensors					
Description	Patient	Patient Weight	Length	Part Number	Quantity
Finger	Adult, Pediatric	≥ 20 kg	1 m	TS-F-D	1
PediTip	Pediatric	15-30 kg		TS-SP-D	1
Fingertip	Adult	> 30 kg		TS-SA-D	
Ear	Adult, Pediatric	≥ 10 kg		TS-E-D	1
Wrap	Adult, Pediatric, Infant	≥ 3 kg		TS-W-D	1
Sensitive Skin	Adult, Pediatric, Infant, Neonate	No limit		TS-SE-3	3

Requires Compatible Interconnect Cable



TS-G3

TruSignal™ Interconnect Cables			
Connector	Length	Part Number	Quantity
GE	3 m	TS-G3	1
Ohmeda		TS-H3	
Datex		TS-N3	



TruSignal™ Disposable Sensors					
Description	Patient	Patient Weight	Length	Part Number	Quantity
Adult/Pediatric	Pediatric, infant	≥ 20 kg	0.3 m	TS-AP-10	10
				TS-AP-25	25
AllFit	Adult, pediatric, infant, neonate	No limit	0.5 m	TS-AF-10	10
				TS-AF-25	25
Adhesive Wrap	Adult	> 20 kg	0.3 m	TS-AAW-10	10
				TS-AAW-25	25
	Pediatric, infant	3-20 kg		TS-PAW-10	10
				TS-PAW-25	25

Requires Compatible Interconnect Cable

TruSignal™ Sensor Accessories

Description	Use With	Size	Part Number	Quantity
Replacement wide adhesive tape	TruSignal™ Sensitive Skin Sensor, TruSignal™ Wrap Sensor	-	OXY-RTW	100
Replacement adhesive tape (blue)	TruSignal™ AllFit Sensor	-	OXY-RT	100
Replacement adhesive tape (teddy bears)			OXY-RTB	
Replacement foam wrap	TruSignal™ Sensitive Skin Sensor TruSignal™ Wrap Sensor	Small	OXY-RWS	24
		Medium	OXY-RWM	
		Large	OXY-RWL	
Infant foam sandal	TruSignal™ Sensitive Skin Sensor	-	OXY-SND	3
Replacement Headband	TruSignal™ Ear Sensor	-	OXY-HB	3



Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Masimo™ SpO2 technology.



LNCS Digit Sensor



DCI-P Digit sensor



LNCS Tip-Clip Ear Sensor

Masimo™ Reusable LNCS and RD SET Sensors

Sensors	Masimo™ Reference	Type	Patient	Patient Weight	Part Number	Quantity
LNCS Digit Sensor	1863	DC-I	Adult	> 30 kg	2027258-001	1
	1864	DCI-P	Pediatric	10-50 kg	2027259-001	
LNCS Tip-Clip Ear Sensor	1895	TC-I	Adult	> 30 kg	2027261-001	1
RD SET SpO2 Sensor	4050	DCI	Adult	> 30 kg	MASIMO-RD-4050	1
	4051	DCI-P	Pediatric	10-50 kg	MASIMO-RD-4051	

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.



2027263-002

Masimo™ LNCS & RD SET Interconnect Cables

Description	Size	Masimo™ Reference	Part Number	Quantity
LNC-10 Interconnect Cable	For GE DINAMAP	1814	2027263-001	1
	For GE	2016	2027263-002	1
RED Interconnect Cable	LNC-04	2055	2059005-001	1
	LNC-10	2056	2059006-001	1
	LNC-14	2057	2059007-001	1
RD SET Interconnect Cable	GE-05 (For GE)	4084	MASIMO-RD-4084	1
	GE-12 (For GE)	4085	MASIMO-RD-4085	1
RD SET Interconnect Cable	MD14-05, 14-pin	4080	MASIMO-RD-4080	1
	MD14-12, 14-pin	4081	MASIMO-RD-4081	1
	MD20-05, 20-pin	4103	MASIMO-RD-4103	1
	MD20-12, 20-pin	4104	MASIMO-RD-4104	1

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries



2051070-001

Masimo™ LNCS & RD SET Disposable Sensors

Description	Patient Weight	Masimo™ Reference	Part Number	Quantity	
LNCS SpO2	Adult	>30 kg	1859	2027253-001	20
	Pediatric	10-50 kg	1860	2027254-001	20
	Infant	3-20 kg	2328	2051068-001	20
	Neonatal PT	< 1 kg	2330	2051070-001	20
	Adult / Neonatal	> 40 kg / < 3 kg	2329	2051069-001	20
RD SET® SpO2	Adult	>30 kg	4000	MASIMO-RD-4000	20
	Pediatric	10-50 kg	4001	MASIMO-RD-4001	20
	Infant	3-20 kg	4002	MASIMO-RD-4002	20
	Adult / Neonatal	> 40 kg / < 3 kg	4003	MASIMO-RD-4003	20
	Premature Baby	< 1 kg	4004	MASIMO-RD-4004	20

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Masimo™ LNCS & RD SET Sensor Accessories

Description	Masimo™ Reference	Part Number	Quantity	
LNCS Replacement Tape	Infant	2307	2051071-001	100
	Neonatal	2308	2051073-001	100
	Neonatal PT	2309	2051074-001	100
RD Set Replacement Tape	Infant	4046	MASIMO-RD-4046	102
	Neonatal	4047	MASIMO-RD-4047	102
RD Set Replacement Wrap	Neonatal PT	4048	MASIMO-RD-4048	10

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



Nellcor™ Interconnect Cables

Description	Connector	Length	Part Number	Quantity
OxiMax™ SpO2 Interconnect Cable	GE	1.2 m	2021406-002	1
		3 m	2021406-001	1
	GE for Corometrics™ only	2.9 m	2025350-001	1
OxiSmart SpO2 Interconnect Cable	GE	2.9 m	2006644-001	1
Interconnect Cables*	DOC-4	1.2 m	2059000-001	1
	DOC-8	2.4 m	2059001-001	1
	DOC-10	3 m	2008773-001	1

* These accessories are sold ONLY with the equipment

414248-001



70124033



70124021



70124035

Nellcor™ Reusable Sensors

Sensors	Patient	Patient Weight	Nellcor™ Reference	Part Number	Quantity
DuraSensor Finger Probe	Adult	< 30 kg	DS100A-1	70124021	1
Multisite SpO2 Sensor	Adult	> 1 kg	D-YS	70124033	1
Two-piece SpO2 sensor	Neonatal / Adult	< 3 kg / > 40 kg	OXI-A/N	70124035	1 (delivered with 24 adhesive strips)
	Pediatric, infant	3-40 kg	OXI-P/I	414248-001	1 (delivered with 24 adhesive strips)

Requires Compatible Interconnect Cable

These accessories are sold ONLY with the equipment

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



70124027



70124032

Nellcor™ Disposable Sensors

Sensors	Nellcor™ Reference	Patient	Patient Weight	Part Number	Quantity
Adhesive SpO2 sensor	MAX-A	Adult	> 30 kg	70124027	24
Long Finger Adhesive sensor	MAXAL			2028117-001	24
Adhesive SpO2 sensor	MAXP	Pediatric	10-50 kg	70124022	24
	MAXI	Infant	3-20 kg	70124026	24
	MAX-N	Adult	> 40 kg	70124032	24
	Neonate	< 3 kg			

Requires Compatible Interconnect Cable
 These accessories are sold ONLY with the equipment



Nellcor™ Fixation Clip

Description	Nellcor™ Reference	For Use With	Part Number	Quantity
SpO2 Ear-Clip	D-YSE	D-YS	70124034	1

These accessories are sold ONLY with the equipment




Nellcor™ Sensor Adhesive Tape

Product	Patient	Nellcor™ Reference	Use With	Part Number	Quantity
Sensor Wrap Replacement Tape	Adult/Neonatal	ADH-A/N	OXI-A/N, D-YS	2016130-001	100
	Pediatric/Infant	ADH-P/I	OXI-P/I, D-YS	2016131-001	100

These accessories are sold ONLY with the equipment

Not available in all countries

Temperature



2104178-001



2104178-002



2104179-001

Temperature Cables				
Description	Connection	Length	Part Number	Quantity
400 Series for Reusable Probes	Single	0.5 m	2104178-001	1
	Dual	0.5 m	2104178-002	1
400 Series for disposable Probes	Single	3.6 m	2104179-001	1

Product compatibility

- CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Not available in all countries



2107176-032



2107176-033

Temperature Probes Reusable					
Description	Connection	Matching Code on the Label	Length	Part Number	Quantity
Esophageal/rectal	Adult	1010132011GE	280 cm	2107176-031	1
	Pediatric	1010132010GE	280 cm	2107176-032	1
Skin	Adult	1010132023GE	280 cm	2107176-033	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 3000/4000/5000 patient monitors
- Solar™ 8000M with TRAM module

Specifications available on request

Not available in all countries



Invasive Blood Pressure



Argon™ - Single



Edwards™ LS TruWave™ - Single



Argon™ - Dual



Edwards™ LS TruWave™ - Dual

Invasive Blood Pressure Cables

Description	Connection	Length	Part Number	Quantity
Argon™	Single	3.6 m	2104166-001	1
		1.2 m	2104166-002	
	Dual	3.6 m	2104166-003	1
		1.2 m	2104166-004	
Edwards™ LS TruWave™	Single	3.6 m	2104162-001	1
		1.2 m	2104162-002	
	Dual	3.6 m	2104162-003	1
		1.2 m	2104162-004	

Product compatibility


- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-SPMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- CARESCAPE ONE Monitor
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

Not available in all countries



INVASIVE BLOOD PRESSURE CABLES



ICU Medical Transpac-IV® - Dual



Utah/Biotrans™ - Single



Utah/Biotrans™ - Dual



Dual Adapter Cable

Invasive Blood Pressure Cables				
Description	Connection	Length	Part Number	Quantity
ICU Medical Transpac-IV®	Single	3.6 m	2104158-001	1
		1.2 m	2104158-002	
	Dual	3.6 m	2104158-003	1
		1.2 m	2104158-004	
Utah/Biotrans™	Single	3.6 m	2104107-001	1
		1.2 m	2104107-002	
	Dual	3.6 m	2104107-003	1
		1.2 m	2104107-004	
Dual Adapter Cable	-	30 cm	2103966-001	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

Not available in all countries



Cuffs






Hoses					
Description	Patient	Compatibility	Length	Part Number	Quantity
 <p>2058203-001 2058203-002 2058203-003</p>	Adult/ Pediatric	V100/VC150, DINAMAP Pro Series, DINAMAP Procure	3.6 m	2058205-001	1
			7.2 m	2058205-002	
		Block style monitor connector for Coro, PDM and some older legacy monitors	1.2 m	2058203-001	
			3.6 m	2058203-002	
			7.2 m	2058203-003	
		CARESCAPE ONE	1.2 m	2087389-001	
			3.6 m	2087389-002	
			7.2 m	2087389-003	
		 <p>2058204-002 2058204-003</p>	Neonate	CARESCAPE ONE	
3.6 m	2089791-003				
Block style monitor connector for Coro, PDM and some older legacy monitors	2.4 m			2058204-002	
	3.6 m			2058204-003	
V100/VC150, DINAMAP Pro Series, DINAMAP Procure	3.6 m			2058206-001	

Not available in all countries

DINACLICK










CONNECTORS

	PART NUMBER	ACCESSORIES	QUANTITY/PACKAGE
	2059294-001	2T DINACLICK FEMALE RECEPT	10/PACK
	2055630-001	1T DINACLICK 1/8" TUBE FEMALE RECEPT	10/PACK
	2055630-002	1T CONN 2-1 DINACLICK 5/32" TUBE FEMALE RECEPT	10/PACK

WALL MANOMETERS AND ACCESSORIES

PART NUMBER	WALL MANOMETER AND ADAPTERS FOR MANUAL INFLATION SYSTEMS	QUANTITY/PACKAGE
2059289-001	INCLUDES COILED TUBING, ADAPTER WITH DINACLICK CONNECTOR, BULB AND VALVE	5/BOX
2059289-002	INCLUDES ADAPTER WITH DINACLICK CONNECTOR AND 13" TUBING (NO COILED TUBING, NO BULB AND VALVE)	5/PACK

ADAPTERS

	PART NUMBER	DESCRIPTION	QUANTITY/PACKAGE
	2059291-001	Y-ADAPTER: 1T MALE BAYONET TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059291-002	Y-ADAPTER: 1T FEMALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059291-003	Y-ADAPTER: 1T MALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059291-004	Y-ADAPTER: 1T UNFITTED TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-001	X-ADAPTER: 2T MALE BAYONET TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-002	X-ADAPTER: 2T FEMALE SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-003	X-ADAPTER: 2T MATED SUBMIN TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059292-004	X-ADAPTER: 2T MALE SCREW TO 2T DINACLICK FEMALE RECEPT	10/PACK
	2059295-001	WASH PLUGS – FOR NEO-SNAP FEMALE & DINACLICK MALE	50/PACK

Not available in all countries





Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It is economically priced to help support your infection control efforts.



20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK	1-TUBE BAYONET CONNECTOR
Infant	Orange	8-13 cm	SFT-P1-2A-5S	2360
Child	Green	12-19 cm	SFT-P2-2A-5S*	2361*
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S	
Small adult	Royal Blue	17-25 cm	SFT-A1-2A-5S*	2362*
Small adult Long	Royal Blue	17-25 cm	SFT-A1-2A-L-5S	
Adult	Navy	23-33 cm	SFT-A2-2A-5S*	2363*
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S*	2140*
Large Adult	Pink	31-40 cm	SFT-A3-2A-5S*	2364*
Large Adult Long	Pink	31-40 cm	SFT-A3-2A-L-5S	2141
Thigh	Brown	38-50 cm	SFT-T1-2A-5S	2365

*Part number with an "-INT" suffix, denoting that the Instructions For Use (IFU) is available in multiple languages.

ASSORTMENTS DINACLICK CONNECTOR			2-TUBE DINACLICK
PRODUCT			
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult			2059303-501
20 Cuffs : 6 Infant, 10 Child, 4 Small Adult			2059303-502
20 Cuffs : 3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm)			2059303-503
20 Cuffs : 2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh			2059303-504
4 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF			2059303-506

Not available in all countries



CLASSIC-CUF

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong enough to withstand repeated inflations – yet economical enough for single patient- use.



5 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK	1-TUBE BAYONET CONNECTOR
Infant	Orange	8-13 cm	CLA-P1-2A-5S	2300
Child	Green	12-19 cm	CLA-P2-2A-5S	2301
Small Adult	Light Blue	17-25 cm	CLA-A1-2A-5S	2302
Adult	Navy	23-33 cm	CLA-A2-2A-5S	2303
Adult long	Navy	23-33 cm	CLA-A2-2A-L-5S	2306
Large Adult	Pink	31-40 cm	CLA-A3-2A-5S	2304
Large Adult Long	Pink	31-40 cm	CLA-A3-2A-L-5S	2307
Thigh	Brown	38-50 cm	CLA-T1-2A-5S	2305



ASSORTMENTS DINACLICK CONNECTOR

PRODUCT

3 Cuffs : 1 Small Adult, 1 Adult, 1 Large Adult

20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult

2-TUBE DINACLICK

2059304-502

2059304-501

Not available in all countries





CRITIKON™ DURA-CUF™ is gentle and designed for durability. Constructed with materials not made with PVC or rubber latex.

Easy-to-use, fast **DINACLICK connector** enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors. The audible click lets users know there is a proper connection.



5 cuffs/box



SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK	1-TUBE BAYONET CONNECTOR
Infant	Orange	8-13 cm	DUR-P1-2A-5S	002274
Child	Green	12-19 cm	DUR-P2-2A-5S	002275
Small Adult	Royal Blue	17-25 cm	DUR-A1-2A-5S	002276
Small Adult Long	Royal Blue	17-25 cm	DUR A1-2A-L-5S	
Adult	Navy	23-33 cm	DUR-A2-2A-5S	002277
Adult long	Navy	23-33 cm	DUR-A2-2A-L-5S	002130
Large Adult	Wine	31-40 cm	DUR-A3-2A-5S	002278
Large Adult Long	Wine	31-40 cm	DUR-A3-2A-L-5S	002131
Thigh	Brown	38-50 cm	DUR-T1-2A-5S	002279

ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
3 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult	2059301-501
6 cuffs : 2 Infant, 3 Child, 1 Small Adult	2059301-502
6 cuffs : 1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult	2059301-503
6 cuffs : 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh	2059301-504

Not available in all countries

YELLOW CLASSIC-CUF

Yellow CRITIKON™ CLASSIC-CUF™ is designed for single patient use within an isolation environment. These vinyl cuffs are priced as disposable-use cuffs and support your infection prevention protocols. They provide the same quality and assurance that you expect from the entire line of CRITIKON cuffs, while tailored for isolation areas.



20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACCLICK
Infant	Yellow/Blue	8-13 cm	CLA-P1-2A-X-5S
Child	Yellow/Blue	12-19 cm	CLA-P2-2A-X-5S
Small Adult	Yellow/Blue	17-25 cm	CLA-A1-2A-X-5S
Adult	Yellow/Blue	23-33 cm	CLA-A2-2A-X-5S
Large Adult	Yellow/Blue	31-40 cm	CLA-A3-2A-X-5S
Thigh	Yellow/Blue	38-50 cm	CLA-T1-2A-X-5S

Not available in all countries





GE Healthcare offers the first BP forearm cuff clinically validated for use on obese or overly muscular adult patients. It provides proper fit and comfort to patients, while giving reliable blood pressure readings that support quality clinical care.

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
OBESE ADULT	PURPLE	26-36 CM (FOREARM)	SFT-F1-2A-5S



Not available in all countries

NEONATAL



Providing big benefits for your tiniest patients, GE Healthcare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl.

The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single-patient-use.

20 cuffs/box


			2-TUBE NEO-SNAP		1-TUBE NEO-SNAP CONNECTOR	
SIZE	COLOR	LIMB CIRCUMFERENCE	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF	SOFT-CUF
1	Orange	3-6 cm	CLA-N1-2B	SFT-N1-2B	CLA-N1-1B	SFT-N1-1B
2	LT. Blue	4-8 cm	CLA-N2-2B	SFT-N2-2B	CLA-N2-1B	SFT-N2-1B
3	Green	6-11 cm	CLA-N3-2B	SFT-N3-2B	CLA-N3-1B	SFT-N3-1B
4	Navy	7-13 cm	CLA-N4-2B	SFT-N4-2B	CLA-N4-1B	SFT-N4-1B
5	Pink	8-15 cm	CLA-N5-2B	SFT-N5-2B	CLA-N5-1B	SFT-N5-1B

NEONATAL ASSORTMENT PACK	CUFF TYPE	PART NUMBER	DESCRIPTION
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-003	1 Neo #3, 1 Neo #4, 1 Neo #5 (3 cuffs/pack)
1-Tube Neo-Snap Connectors	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)


Not available in all countries



Cardiac Output



CARDIAC OUTPUT



Cardiac Output Cables

Description	Length	For use with	Part Number	Quantity
Cardiac Output Cable	3.6 m	CARESCAPE™ Monitors B450, B650 and B850 with the CARESCAPE™ Patient Data Module (PDM)	2104180-001	1
	2.4 m		2104180-002	1

Specifications available on request

Not available in all countries

OTHER ACCESSORIES

Other Accessories				
Description	For use with	Manufacturer Reference	Part Number	Quantity
CO Set for Room Temperature Injectate, disposable	E-COP, E-COPsv	93610	16574VS	1
Catheter Connecting Cable	E-COP, E-COPsv	COM2-CC	16590-HEL	1
Thermodilution Catheter, disposable	E-COP, E-COPsv	131HF7	M1238258	5

These accessories are sold ONLY with the equipment

PiCCO™ CABLES

PiCCO™ Cardiac Output Cables				
Description	Length	For use with	Part Number	Quantity
Continuous Cardiac Output Cable	4.1m	E-PiCCO module	M1154585	1
Injectate Sensor Cable	0.9m		M1182448	1

These accessories are sold ONLY with the equipment



GE Healthcare compatible Cardiac Output probes are available from multiple manufacturers. You can purchase the following product portfolios from our partners directly:

- Vyair Medical
- Bath Probes
- Argon Medical
- Spectramed Temperature Probes
- Edwards Life Sciences
- Injectate Bath Probes
- Flow-through Injectate Temperature Probes
- Getinge
- PiCCO Accessories

Not available in all countries

Adequacy of Anesthesia



M1038681

M1174413



ENTROPY

Sensors			
Product	Use	Part Number	Quantity
Entropy™ EasyFit Sensor	Disposable	M1174413	25
Entropy™ Sensor		M1038681	

Cable			
Product	Length	Part Number	Quantity
Entropy™ Cable	3.5 m	M1050784	1

Product compatibility

E-ENTROPY module

The GE Healthcare Entropy™ module, E-ENTROPY, and accessories are indicated for adult and pediatric patients older than 2 years within a hospital for monitoring the state of the brain by data acquisition of electroencephalograph (EEG) and frontal electromyograph (FEMG) signals during general anesthesia.

Not available in all countries



Electromyography (EMG)
Requires 5 electrodes



Kinemyography (KMG)
Requires 2 electrodes



888418/897439



888416

Sensors				
Product	Use	Patient	Part Number	Quantity
NMT MechanoSensor	Reusable	Adult	888418	1
		Pediatric	897439	
NMT ElectroSensor		Adult, Pediatric	888416	1
NMT Regional Block Adapter		Adult, Pediatric	888417	1

Accessories				
Product	Use	Patient	Part Number	Quantity
NMT Electrodes	Disposable	Adult, Pediatric	57268-HEL	30

Cable			
Product	Length	Part Number	Quantity
NMT Sensor Cable	3 m	888414	1


Product compatibility

E-NMT module

Not available in all countries



Neurology



898050

Reusable Preconfigured EEG Leadsets

Product	Length	Part Number	Quantity
Basic montage	0.6 m	898050	1
AEP montage		898052	1

Not available in all countries



EEG ELECTRODES



2083319-001



2083320-001

Disposable EEG Electrodes

Product	Part Number	Quantity
EEG Clip electrodes Ag/AgCl	M1240794	40 pouches of 25 (1,000 total)
	572686	500
EEG Prewired cup electrodes	2083319-001	46 pouches of 10 (460 total)
EEG Prewired wet gel electrodes	2083320-001	42 pouches of 10 (420 total)

All electrodes within a set-up should be similar

Product compatibility

E-EEG(X) module and N-EEG(X) headbox

AEP Accessories

Product	Part Number	Quantity
Disposable ear tips for use with earphone set	8001233	20

To use with 898052

Product compatibility

E-EEG(X) module and N-EEG(X) headbox

Not available in all countries



Patient Spirometry



896952



8001948



733910-HEL



73393

Sensors				
Description	Type	Approved with Modules	Part Number	Quantity
D-lite™ +	Single Patient Disposable	E-COV, E-COVX, E-CAiOV, E-CAiOVX, E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX	896952	50
Pedi-lite+			8001948	
D-lite™	Re-usable	E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX	733910-HEL	1
Pedi-lite			73393	

Specifications available on request

(1) Not available in the US

Not available in all countries



SPIROMETRY



890031

Spirometry Tubes

Description	Length	Type	Approved with Modules	Part Number	Quantity
Spirometry Tubes	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	890031	5
	3m			884101	
	6m			891191-HEL	



D-lite™ ++ set

Kits

(include gas sample line, sensor and spirometry tube)

Description	Length	Type	Approved with Modules	Part Number	Quantity
D-lite™ ++ set	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX	2104297-001	20
	3m			2104297-002	
D-lite™ + kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8004381	50
	3m			8004382	
Pedi-lite+ kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8002718	50
	3m			M1032634	

Specifications available on request

Pedi-lite+ kit

(1) Not available in the US

Not available in all countries



Gas Monitoring



ICU AIRWAY ADAPTER SETS



Airway Gas Sampling Sets					
Description	Type	Approved with Modules	Length	Part Number	Quantity
Adult/Pediatric	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX	3 m	8090121531	20 pcs/box
			2 m	8090121521	
Infant		E-sCO, E-sCOV, E-sCOVX	2 m	8090121522	20 pcs/box

Not available in all countries

GAS SAMPLING LINES AND ADAPTERS



Gas Sampling Lines					
Description	Type	Approved with Modules	Length	Part Number	Quantity
PVC/PE gas sampling line	Single Patient Disposable	E-miniC, E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOVX, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹ , N-CAiO	2 m	2097307-001	10
			3 m	2097307-002	
			6 m	2097307-003	5



3090121002

Gas Monitoring Airway Adapter				
Description	Type	Approved with Modules	Part Number	Quantity
Airway Adapter with Luer connector	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹ , N-CAiO	3090121002	25

WATER TRAPS



881319-HEL



M1200227



876446-HEL



M1182629



8002174

Water Traps – ICU				
Description	Type	Approved with Modules	Part Number	Quantity
D-fend™+	Single Patient Disposable	E-CO, E-COV, E-COVX	881319-HEL	10
D-fend™ Pro+		E-sCO, E-sCOV, E-sCOVX	M1200227	
Mini D-fend™		E-miniC	8002174	

Specifications available on request


Water Traps – Anesthesia				
Description	Type	Approved with Modules	Part Number	Quantity
D-fend™	Multi-Patient	E-CAiO, E-CAiOV, E-CAiOVX	876446-HEL	10
D-fend™ Pro		E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹	M1182629	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

Not available in all countries



ANESTHESIA & RESPIRATORY CARE



ANESTHESIA & RESPIRATORY CARE

Absorbents

- Absorbents
- Canisters
- Dust Filter

Vaporizer Bottle Adapters

Filters & Breathing Masks

Reservoir Bags & Gas Scavenging

- Anesthesia Breathing Bags
- Scavenging Bags, Tubes, & Connectors

Patient Breathing Circuits

- Circuits
- Circuit Components
- Circuit Kits
- Circuit Startup Kits

Flow Sensors - Respiratory

- Cables
- Heater
- Respiratory Flow Sensors and Expiratory Valve Assemblies

Flow Sensors - Anesthesia

Nebulizer Accessories


Anesthesia O2 Cells

Respiratory Safety Guard

Pressure Catheters



Anesthesia & Respiratory care



Disposable CO2 Absorbers			
Description	Compatibility	Part Number	Quantity
AMSORB® Plus G-Can Absorber, 1.0L	All Aisys™, Avance™, Aespire™, 9100c and ADU variants	2105489-003	8
AMSORB® Plus Prefilled Cartridge, 1.2L	All Aestiva™ variants	2105489-004	12
AMSORB® Plus CareCan Absorber, 1.4L	Carestation™ 600 series & 750 series*	2105489-006	8

Reusable CO2 Absorbent Components			
Description	Compatibility	Part Number	Quantity
AMSORB® Plus Jerican Loosefill Absorbent, 5.0L	GE Healthcare Reusable Canisters	2105489-007	2
ABS Reusable CO ₂ Canister	All Aisys™, Avance™, Aespire™ & 9100c variants	1407-7004-000	1
CBS Reusable CO ₂ Canister	Carestation™ 600 series & 750 series*	2071165-001-S	1
Foam Dust Filter	ABS Reusable Canister	2105489-001	40

Easy-Fil Vaporizer Bottle Adapters			
Vaporizer Compatibility	Color Indicator	Part Number	Quantity
Isoflurane Easy-Fil	Purple	2106572-001	1
Halothane Easy-Fil	Red	2106572-002	
Sevoflurane Easy-Fil	Yellow	2106572-003	



*Not cleared approved by the US FDA. Not for sale in the United States.

FILTERS AND BREATHING MASKS



HMEFs and Filters – Disposable				
Description	Size	Sampling Port	Part Number	Quantity
HMEF	Neonatal	Included	2106570-006	50
	Pediatric		2106570-009	
	Adult		2106570-010	
Bacterial and Viral Filter	Adult	Included	2106570-007	50
HEPA Filter	Adult		2106570-008	



Anesthesia Breathing Masks – Disposable			
Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-012	50
	1 – Small Child	2106570-013	
22mm	2 – Child	2106570-014	50
	3 – Large Child	2106570-015	
	4 – Adult	2106570-016	50
	5 – Large Adult	2106570-017	
6 – Extra Large Adult	2106570-018		

Anesthesia Breathing Masks – Reusable			
Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-011	1
		2106570-029	10
	1 – Small Child	2106570-019	1
		2106570-030	10
22mm	2 – Child	2106570-020	1
		2106570-031	10
	3 – Large Child	2106570-021	1
		2106570-032	10
	4 – Adult	2106570-022	1
		2106570-033	10
	5 – Large Adult	2106570-023	1
		2106570-034	10

Not available in all countries



RESERVOIR BAGS AND GAS SCAVENGING

Anesthesia Manual Breathing Bags – Reusable



Bag Connection Size	Material	Bag Size	Part Number	Quantity
22mm	Clear Silicone	0.5L	2106570-001	1
			2106570-025	5
		1.0 L	2106570-002	1
			2106570-026	5
		2.0 L	2106570-003	1
			2106570-027	5
		3.0 L	2106570-004	1
			2106570-028	5

Passive Scavenging Tubes – Disposable



Connection	Length	Type	Part Number	Quantity
30M-30M	5 m	Tube	2105488-006	1
			2105488-022	10

Gas Scavenging Bag – Reusable



Size	Compatibility	Connection	Part Number	Quantity
3.0L	Aisys, Avance, Aespire™ (ABS Breathing Systems)	30M	2106570-024	1

Straight Connectors – Scavenging

Connection	Reprocessable	Part Number	Quantity
30M-19M	Disposable	M1003947	1
		1500-3376-000	5
30M-309M		M1003134	1

Not available in all countries



PATIENT BREATHING CIRCUIT KITS – START UP

All startup kits include patient breathing circuit, reservoir bag tube, reservoir bag, gas sampling lines and gas sampling exhaust lines. These are primarily intended for initial anesthesia machine setup.



Patient Breathing Circuit Kits – Start Up – Adult

Type	Gas Sampling Line	Gas Sampling Exhaust Line	Limb Circuit	Breathing Circuit	Breathing Bag		Part Number	Quantity
				Length	Manual Tube 1.5 M	Size		
Disposable	3m PVC/PE	0.18m & 1m	Dual	Expandable up to 3 m	Straight	3L	2096533-001	1
						2L	2096533-003	1
			Single Split Lumen	1.8 m	Angled*	2L	2096536-002	2
					Straight	2L	2096536-003	2
					Angled*	2L	2105488-019	1
Reusable			Dual	1.5 m	Straight	2L	2096534-009	1
					Angled*	2L	2096536-001	1

*For use with bag support arm

Not available in all countries



PATIENT BREATHING CIRCUITS

Our wide range of patient breathing circuits and components offer robust and high quality options for both disposable and reusable applications.

Patient Breathing Circuit Kits

Type	Size	3m PVC/PE Gas Sampling Line	Circuit Type	Breathing Circuit	Breathing Bag		Includes	Part Number	Qty
				Length	Manual Tube	Size			
Disposable	Adult	Yes	Dual Limb	Expandable up to 3 m	Not included	3L	-	2096535-006	20
				1.8 m	Angled*, 1 M		-	2105488-013	20
					Straight, 1.5 M		-	2105488-014	20
		No	Single Split Lumen	1.8 m	Straight	3L	-	2105488-010	20
				Dual Limb	Expandable up to 1.8 m	Not included	-	2105488-011	20
			1.5 m		Not included	-	2x Water Traps	2096535-001	20
	Straight, 0.8 M			-	1x Water Trap	2105488-017	20		
	Neonatal	Yes	Dual Limb	1.5 m	Not included	1L	HME & Filter	2105488-001	20
		No						2105488-021	1
	Reusable	Adult	No	Dual Limb	1.5 m	Not included	2x Water Traps	2096534-005	1
Pediatric		2096534-006						1	

*For use with bag support arm



Not available in all countries

PATIENT BREATHING CIRCUITS



2096534-007
Reusable, Pediatric



2105488-020
Disposable, Pediatric




2105488-015
Disposable, Adult



2105488-009
Disposable, Adult, Inteliflo® Single Split Lumen

Patient Breathing Circuits

Type	Size	Circuit Type	Breathing Circuit Length	Patient Elbow	Gas Sampling Port Location	Part Number	Quantity
Reusable	Adult	Dual Limb	1.5 m	Double Swivel, Removable	Elbow	2096534-007	1
	Pediatric					2096534-008	1
Disposable	Adult	Single Split Lumen	1.8 m	Fixed, Removable	Elbow	2105488-009	20
		Dual Limb	1.5 m	Double Swivel, Removable	Patient Wye	2096535-002	20
	Fixed, Removable			Elbow	2105488-015	20	
	Pediatric	Dual Limb	1.5 m	Double Swivel, Removable	Patient Wye	2105488-020	20



Non-Recirculating Breathing Circuits – Jackson Rees

Type	Breathing Circuit Length	Bag Size	Gas Sampling Port Location	Part Number	Quantity
Disposable	0.8 m	1 L	Elbow	2105488-008	20

Not available in all countries



Reusable Patient Breathing Circuit Tubes

Size	Material	Length	Connections	Part Number	Quantity
Adult	Hytrel®	0.35 m	22F-22F (90 Deg)	2096534-001	1
		0.75 m		2096534-002	1
		0.9 m		2096534-003	1
		1.5 m		2096534-004	1
	Silicone	1 m	22F-22F	2096534-101	1
		1.5 m		2096534-102	1
Pediatric	Silicone	1.5 m	22F-10F	2096534-103	1
	Hytrel®	1.1 m	22F-15F	2105488-007	1
				2105488-023	10



Not available in all countries



PATIENT BREATHING CIRCUITS



Patient Breathing Circuit Components				
Description	Size	Reprocessable	Part Number	Quantity
Y-Piece Connector Kit, 15M-10F	Neonatal	Reusable	2096534-104	10
Patient Angled Y-Piece			2105488-002	10
Patient Straight Y-Piece - Swivel			2105488-004	10
Straight Connector - CBS Bag Arm, 22M-22M	Adult	Reusable	2096533-004	1
ABS Bag Arm Connector, 22M-22M	All		2105488-005	10
Patient Elbow - Swivel	Adult	Reusable	2105488-003	10
Patient Elbow		Disposable	2105488-018	50
Gas Return Adapter, 22M-22F	Adult	Disposable	2105488-016	10

Not available in all countries





Respiratory Flow Sensors and Expiratory Valve Assemblies (EVA)				
Description	Reprocessable	Patient Type	Part Number	Quantity
Flow Sensor Only	Reusable	All	1505-3231-000	1
Whole Assembly (housing and flow sensor)			1505-3848-000	1
Housing Only (includes water trap)			1505-8568-000	1
Flow Sensor Only (Neonatal)		Neonatal	1505-3272-000	1
Flow Sensor Only	Single Patient Use	All	1505-3231-001	1
Whole Assembly (housing and flow sensor)			1505-3848-001	1
Housing Only (includes water trap)			1505-8568-001	1
Flow Sensor Only			1505-3231-010	10
Whole Assembly (housing and flow sensor)			1505-3848-010	10
Housing Only (includes water trap)			1505-8568-010	10

Product compatibility

CARESCAPE™ R860 and Engström Ventilators

Not available in all countries

FLOW SENSORS – RESPIRATORY



M1200693



M1188723

Respiratory Flow Sensor and Expiratory Valve Components

Description	Patient Type	Part Number	Quantity
Flow Sensor Cable	Neonatal	1505-5604-000	1
Exhalation Valve Heater	All	M1200693	1
Exhalation Valve Heater Cable		M1188723	1


Product compatibility

CARESCAPE™ R860 and Engström Ventilators

FLOW SENSORS – ANESTHESIA



2087640-001-S



2089610-001-S



2069358-001-S

Anesthesia Flow Sensors

Machine Compatibility	Autoclavable	Part Number	Quantity
Aisys™, Avance™, Aespire™, Aestiva™	Yes	2087640-001-S	1
	No	2089610-001-S	
Carestation™ 600 series & 750 series*	Yes	2069358-001-S	1

*Not cleared approved by the US FDA. Not for sale in the United States.

Not available in all countries



NEBULIZER ACCESSORIES



Aerogen® Nebulizer Accessories			
Description	Patient Type	Part Number	Quantity
Aerogen® Pro Nebulizer kit with adult and ped tee and cable	Adult and Pediatric	1505-3846-000	1
Nebulizer head with filler cap	All	AG-AP1000	
Nebulizer Cable		1505-5602-000	
Nebulizer tee with plug	Adult	AG-AP1010	5
	Pediatric	AG-AP1020	
	Neonatal	AG-AP1035	
Nebulizer replacement kit		AG-AP1100	1
Aerogen® Solo Nebulizer	Adult	AG-AS3100-S	5
		AG-AS3200-S	10
Nebulizer tee	Adult	AG-AS3010	10
	Pediatric	AG-AS3020	
	Neonatal	AG-AS3035	

Product compatibility

CARESCAPE™ R860 and Engström Ventilators

These accessories are sold ONLY with the equipment.

ANESTHESIA O2 CELLS

Anesthesia O2 Cells			
Description	Machine Compatibility	Part Number	Quantity
O2 Sensor ONLY	Aisys, Avance, Aespire™	M-10	1
O2 Sensor and cable		2082006-001	
O2 Sensor with Adapter and O-Ring	Carestation™ 600 series & 750 serie*	2071018-001-S	



2071018-001-S

*Not cleared approved by the US FDA. Not for sale in the United States.

Not available in all countries



RESPIRATORY INSPIRATORY SAFETY GUARD



Respiratory Inspiratory Safety Guard			
Machine Compatability	Patient Type	Part Number	Quantity
CARESCAPE™ R860	All	5830446	1
		5514015	10

PRESSURE CATHETERS



Respiratory Pressure Catheters			
Description	Machine Compatability	Part Number	Quantity
Intratracheal Pressure Catheter	CARESCAPE™ R860 and Engström Ventilators	M1045564	10

Not available in all countries



DIAGNOSTIC CARDIOLOGY



DIAGNOSTIC CARDIOLOGY

MAC™, CASE™ and CardioSoft™

ECG

- Accessories
- Cables
- Electrodes
- Leadwires
- Recording Paper

Holter Monitors

ECG – CM 4000

- Leadwire Sets
- Accessories

ECG – SEER® 12

- Leadwire Sets
- Accessories

ECG – SEER® 1000

- Leadwire Sets
- Accessories

ECG – SEER® Light

- Leadwire and Cable Kits
- Accessories

MAC™, CASE™ and CardioSoft™



CABLES AND LEADWIRES – IEC

ECG Trunk Cables – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™	
			VU360	7	5500	3500	2000	800			600
Trunk Cable 10-lead IEC	220 cm	2104722-001		X		X	X	X	X		



Integrated Cable and Leadwire Set – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™	
			VU360	7	5500	3500	2000	800			600
ECG Value Cable 10-lead Cable / Leadwire IEC	Chest leads: 80 cm Limb leads: 125 cm Trunk cable: 220cm	2104726-001		X		X	X	X	X		




Not available in all countries



LEADWIRES – IEC

ECG Leadwire Set – Banana – IEC

Description	Length	Part Number	MAC™							CASE™	CardioSoft™
			VU360	7	5500	3500	2000	800	600		
Leadwire Set, Base 10, banana, IEC	Limb leads: 125 cm	2104724-002	X	X	X	X	X	X	X	X	X
Leadwire Set, 4 Add-on, banana, IEC-HEIM	Chest leads: 80 cm	2104750-001	X		X					X	X



ECG Leadwire Set – Grabber – IEC

Description	Length	Part Number	MAC™							CASE™	CardioSoft™
			VU360	7	5500	3500	2000	800	600		
Base 10, IEC	Limb lead: 125 cm	2104748-001	X	X	X	X	X	X	X	X	X
4 Add-on, IEC-HEIM	Chest lead: 80 cm	2104725-002	X		X					X	X



Not available in all countries



CABLES AND LEADWIRES – IEC



2104729-001



2104767-001

Accessories		MAC™								CASE™	CardioSoft™	Quantity
Description	Part Number	VU360	7	5500	3500	2000	800	600				
		Adapter Clip, Wide Mouth Clip	2104729-001	X	X	X	X	X				X
ECG Leadwire Port Plug	2104767-001	X		X					X	X	4	

CARDIOSOFT™



CardioSoft™ Accessories		
Description	Part Number	Quantity
Supply BP Bag Belt TONOPORT	2104824-001	1

Not available in all countries

ECG RECORDING PAPER



DCAR ECG Recording Thermal Paper

Folding	White patient data area	Paper size	Red grid wide	Queue	MAC™							CASE™	Part Number	Quantity
					VU360	7	5500	3500	2000	800	600			
z-fold	No	A4 *	195 mm	Hole			X	X				X	2104769-001	300 sheets/ pack, 8 packs
			205 mm	Block	X	X			X				2104772-001	150 sheets/ pack, 10 packs
		21.5 cm x 28 cm	200 mm	Hole			X	X				X	2104768-001	300 sheets/ pack, 8 packs
		80 mm x 90 mm	75 mm	Block							X		2104773-001	280 sheets/ pack, 10 packs
		110 mm x 140 mm	100 mm	Block						X			2104775-001	200 sheets/ pack, 10 packs
		Yes	21.5 cm x 28 cm	155 mm	Hole	X		X	X					2104770-001
	Hole				X	X	X	X					2104771-001	150 sheets/ pack, 16 packs

For ECG paper for Patient monitoring, please go to page 6.



REUSABLE ECG ELECTRODES

GE Healthcare reusable bulb electrodes with silver-silver chloride (Ag/AgCl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.



Reusable Bulb Electrode		
Description	Part Number	Case Quantity
Reusable Bulb Electrode	2104783-001	6



Bulb Electrode Accessories			
Description	Compatibility	Part Number	Case Quantity
Clamp Limb AHA Value	MAC™ 7, MAC™ 600, MAC™ 800, MAC™ 2000 ECG System	2104784-001	4/set
Clamp Limb IEC Value		2104785-001	

Specifications available on request

Not available in all countries



Holter Monitors



CM 4000 LEADWIRE SETS



2105803-005

CM 4000 Leadwire Set - IEC

Description	Length	Part Number	Quantity
3-leads	80 cm	2105803-001	1
5-leads		2105803-003	1
7-leads		2105803-005	1

Specifications available on request

Not available in all countries



CM 4000 Accessories		
Description	Part Number	Case Quantity
Carry Pouch Reusable	2054933-014	1
Carry Bag Disposable	2054933-015	1
USB Cable	2054933-013	1

Specifications available on request



2104808-001

SEER® 12 – IEC				
Description	Length	Type	Part Number	Quantity
12-Channel Holter ECG Cable	130 cm	Round	2104808-001	1
	105 cm	Flat	2104809-001	

SEER® 12 ECG cables are compatible with SEER® 12 monitors and CardioMem™ CM 3000-12 BT monitors




2104803-001

SEER® 12 Holter Monitor Accessories		
Description	Part Number	Quantity
SEER® 12 lead portable pouch, waist belt, shoulder strap, 1 battery compartment	2104803-001	1
SEER® 12 Carrying pouch, 2 battery compartment	2014606-297	1

2104803-001 is compatible with SEER® 12, CM 3000, CM 3000-12
 2014606-297 must be purchased with service
 Specification available on request

Not available in all countries



2104811-001

SEER® 1000 Leadwire Set – IEC

Description	Length	Part Number	Quantity
5 leadwire, 2 CH	75 cm	2104810-001	1
7 leadwire, 3 CH		2104811-001	
5 leadwire, 2 CH	105 cm	2104812-001	
7 leadwire, 3 CH		2104813-001	
3 leadwire, 3 CH	75 cm	2104818-001	
	105 cm	2104819-001	

Product compatibility

SEER® 1000 ECG cables are compatible with SEER® 1000 24-hour Holter monitor, SEER® 1000 48-hour Holter monitor and SEER® 1000 7-day Holter monitors

Specifications available on request

Not available in all countries



2067634-073

SEER® 1000 Holter Monitor Accessories

Description	Part Number	Quantity
Pouch with waist belt and shoulder strap	2067634-073	1
Pouch Belt Clip	2067634-166	1

Specification available on request



2067634-166

Not available in all countries

SEER® Light Leadwire Set Kits – IEC

Description	Length	Part Number	Quantity
Carry Pouch and 2 Channel Patient cable	140 cm	2042686-003	1
Carry Pouch and 3 Channel Patient cable		2042686-004	
Carry Pouch and 2 Channel Patient cable	60 cm	2042686-007	1
Carry Pouch and 3 Channel Patient cable		2042686-008	

Product compatibility

SEER® Light and SEER® Light Extend

Specifications available on request

Not available in all countries



SEER® Light Accessories		
Description	Part Number	Quantity
Compact Flash	2008750-003	1
Connect Kit	2020044-292	1
Controller Kit	2042685-007	1

Product compatibility

SEER® Light and SEER® Light Extend
Specifications available on request



SEER® Light Holter Monitor Accessories		
Description	Part Number	Quantity
SEER® Light, shoulder straps, waist belt, portable pouch	2104827-001	1

Product compatibility

2104827-001 is compatible with SEER® Light and SEER® Light Extend
Specification available on request

Not available in all countries



MATERNAL INFANT CARE



MATERNAL INFANT CARE

Maternal Fetal Care

Mini Telemetry

- Recording Thermal Paper
- Remote Event Marker

Novii™

- Cables
- Novii™ Patch
- Tape

Corometrics™ and Mini Telemetry

- Cables
- Transducer Belts
- Transducer Button and Loop

Internal Monitoring

- Cables
- Electrode
- Pad

Infant Care

Giraffe™/Panda™

- Air filters
- Clip
- Cover

Phototherapy accessories

- Pad Covers

Resuscitation

- T Piece

Lullaby™ Warmer


- Temperature Probes
- Temperature Probe Patches

Shuttle™

- Patient Restraint

Maternal Fetal Care

- Mini Telemetry
- Novii™
- Corometrics™



MIC Fetal Monitoring Recording Thermal Paper

Scaling	Corometrics™ 170 series	Corometrics™ 250cx series	Part Number	Quantity
30-240 bpm scaling	X	X	2104907-001	40 packs
50-210 bpm scaling	X	X	2104908-001	

For ECG paper for Patient monitoring and Diagnostic cardiology, go to pages 6 and 58.

Specifications available on request



2107349-001

Accessories		
Description	Part Number	Quantity
Remote Event Marker	2107349-001	1

Not available in all countries





Novii™ Accessories			
Description		Part Number	Quantity
Novii™ Patch		107-PT-004-10	10
		107-PT-004-50	50
Novii™ Interface Cable	GE FECG	105-PT-102	1
	GE MECG	105-PT-104	
	GE UA	105-PT-106	
3M™ Red Dot™ 2236 Skin Prep Tape		100-PT-007	1
FECG/MECG Y Adapter Cable, 120 Series		1442AAO	1


Not available in all countries





MECG Accessories - Vyaire			
Description		Part Number	Quantity
Gray MECG Cable	3 Lead Set, Model 118	1554AAO	1
	3 Lead Set IEC, Model 118	1554BAO	1
120 Green MECG Cable	3 Lead Set	1553AAO	1
	3 Lead Set IEC	1553BAO	1

Sold with equipment only.



2015827-001

Transducer Belts			
Description		Part Number	Quantity
Single patient use	Button-style, blue	2015827-001	10




2108347-001

Ultrasound Transducers		
Description	Part Number	Quantity
Ultrasound Transducer Button and Loop	2108346-001	1
TOCO Transducer Button and Loop	2108347-001	1

Not available in all countries





2109002-001



2109006-001

Fetal Spiral Electrode Accessories

Description	Type	Part Number	Quantity
Cable	Rectangular Connector	2109005-001	1
	Round Connector	2109004-001	1
ATTACHMENT PAD	Large, Single Patient Use	2109006-001	50
SPIRAL ELECTRODE	Single Helix, Single Patient Use	2109002-001	50

For ECG Electrodes for Patient Monitoring, go to pages 14 and 15.

Not available in all countries



Infant Care

- Giraffe™/Panda™
- Phototherapy accessories
- Resuscitation



GIRAFFE™/PANDA™



2107356-001



2107357-001

Accessories		
Description	Part Number	Quantity
Air Filters	2107356-001	10
Bag and mask holder clip	2107357-001	1

Not available in all countries



Giraffe™ Hood Covers

Description	Part Number	Quantity
OmniBed™ Cover	6600-0825-800	1
Incubator Cover	6600-0846-800	1

Sold with equipment only

Not available in all countries





Flat cushioned cover with straps
The flat, cushioned BiliSoft™ Cover lets you swaddle baby and phototherapy pad together.



Cushioned cover with developmental positioning foot roll and straps
The BiliSoft™ Nest offers sick babies the boundaries and support they need.

BiliSoft™ Disposable Pad Covers					
Description	BiliSoft™ 1.0	BiliSoft™ 2.0	Size	Part Number	Quantity
Disposable Pad Covers, Flat		X	Small	2105422-001	20
				2105418-001	50
	X		Small	2107201-001	20
				2107200-001	50
Disposable Nest Pad Covers		X	Small	2105420-001	15
	X			2107202-001	
	X	X	Large	2105421-001	

Product compatibility

The BiliSoft™ Pad Covers and Nest are designed for use with the Bilisoft™ 1.0 and 2.0 Phototherapy System, which is used to treat both premature and full-term infants for hyperbilirubinemia.

Available in small and large sizes to match the size of the Light Pad*

* Note: the BiliSoft™ 2.0 Disposable Pad Cover (Flat and Nest) in Size Large also fits the BiliSoft™ Phototherapy System (Legacy) large pad.

Specifications available on request



Not available in all countries

T-Piece Circuit: Disposable resuscitator designed for neonates, provides accurate reliable ventilation and pulmonary resuscitation.



2109061-001
Disposable T Piece
with Mask

Infant Resuscitator T-Piece Circuit Kit

Description	Mask Size	Part Number	Quantity
Disposable T-Piece	-	2109059-001	20
Disposable T-Piece with Mask	0	2109060-001	10
	1	2109061-001	10

Product compatibility

With a universal connector, it fits with all commonly used manual resuscitator units and does not require an adapter.

Specifications available on request



Panda™ Warmer



StAr

Not available in all countries



LULLABY™ WARMER

Temperature Probes - Vyaire		
Description	Part Number	Quantity
Disposable Temperature Probe	6600-0196-700	50
	6600-0208-700	10
Reusable Temperature probe – Lullaby™ Warmer	2044384-001	1

Sold with equipment only in most countries.

Temperature Probe Patches - Vyaire		
Description	Part Number	Quantity
Patch Heat Reflecting ROHS	0203-1980-300	50
Heat Reflecting Patches	2044385-001	50

Sold with equipment only in most countries.

SHUTTLE™

Transport		
Description	Part Number	Quantity
Patient restraint	M1169540	5

Not available in all countries





Product may not be available in all countries and regions.
Full product technical specifications are available upon request.
Contact a GE Healthcare representative for more information.

Please visit: www.gehealthcare.com/promotional-locations

Data subject to change.

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges.

From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

© 2021 General Electric Company – All rights reserved.

GE, the GE monogram, MAC[™], CASE[™], CardioSoft[™], CardioMem[™], SEER[®], Aespire[™], Aestiva[™], Aisys[™], Carestation[™], CARESCAPE[™], Solar[™], Dash[™], Entropy[™], D-lite[™], D-fend[™], Corometrics[™], Shuttle[™], Panda[™], Lullaby[™] and TruSignal[™] are trademarks of General Electric Company. ICU Medical Transpac is a registered trademark of ICU Medical, Inc.; Utah/Biotrans is a registered trademark of Utah Medical Products, Inc. Edwards Lifesciences TruWave is a trademark of Edwards Lifesciences Corporation. Argon is a trademark of Argon Medical Devices. AMSORB and G-Can are trademarks of Armstrong Medical. Hytrel is a registered trademark of E. I. du Pont de Nemours and Company or its affiliates. Inteliflo is a registered trademark of Flexicare Medical Ltd. Nellcor and Oximax are trademarks of a Medtronic Company. Masimo and SET are trademarks of Masimo Corporation. All other trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

JB01042XE Rev 2