



SpO₂ Connection Guide



Contents

Sensor Compatibility	Page
TruSignal™*	3
Masimo™ RD	7
Masimo LNCS®	8
Nellcor™ OxiMax™	10
Nellcor OxiSmart®	12
Telemetry	14

			TruSignal	Masimo RD	Masimo LNCS	Nellcor Oximax	Nellcor OxiSmart	Telemetry
Pulse oximeters	3775, 3700 Series, 3800, 3900, 3900P		4					
	GE Ohmeda TruSat™* oximeter		3					
	GE Ohmeda TuffSat™** oximeter		4					
Configured patient monitors	Cardiacap™*/5		5					
	Cardiacap/5 with N-XOSAT		4					
	Cardiacap/5 with N-XNSAT						11	
	CARESCAPE™* V100		6	7	8	10		
	CARESCAPE VC150 Monitor		6	7	8	10		
	Dash 2500				8	10		
	Dash 3000/4000/5000			7	8	10	11	
	GE DINAMAP™** ProCare				9	10		
	GE DINAMAP PRO 1000 v3				9	10		
	GE DINAMAP PRO Series v2				9	10		
PROCARE monitors	S/5 Light Monitor		5					
	PROCARE* Monitor B40(i) Monitor		6	7	8	10		
	PROCARE Monitor B20(i) Monitor		6	7	8	10		
Modular patient monitors	CARESCAPE Monitor B850	with CARESCAPE Patient Data Module		7	8	10		
		with E-MASIMO		7	8	10		
		with E-NSAT, E-NSATX						
		with E-PSM(P)	6					
		with E-PRESTN, E-RESTN	5, 6		8	10		
		with Tram module ¹			8		11	
		with CARESCAPE ONE	6		8	10		
	CARESCAPE Monitor B650	with CARESCAPE Patient Data Module		7	8	10		
		with E-MASIMO		7	8			
		with E-NSAT, E-NSATX				10		
	with E-PSM(P)	6						
	with E-PRESTN, E-RESTN	6		8	10			
	with CARESCAPE ONE	6						
CARESCAPE Monitor B450	with CARESCAPE Patient Data Module		7	8	10			
	with E-MASIMO		7	8				
	with E-NSAT, E-NSATX				10			
	with E-PSM(P)	6						
	with CARESCAPE ONE	6		8	10			
B105, B125, B30 Monitor, FM, FM Light Monitors			6	7				
S/5, AS/3, CS/3 modular monitors	with E-MASIMO				8			
	with E-NSAT, E-NSATX				8			
	with E-PSM(P)	6						
	with E-PRESTN, E-RESTN	6						
	with M-NSAT	4					11	
	with M-OSAT	5						
	with hemodynamic M-modules							
Solar™* 8000/8000M/8000i and 9500	with CARESCAPE Patient Data Module			7	8	10		
	with Tram module					10	11	
Other monitoring devices	ApexPro™** telemetry transmitter							12
	CardioLab™**/Mac-Lab* with Tram module					10	11	
	Corometrics™** 118						11	
	Corometrics 120/2120 Series						11	
	Corometrics 250 Series			7	8			
	Responder* 3000						11	

¹ Not compatible with software version 3

* For monitors or modules no longer for sale, applicable sensors are available as service items.

Data subject to change.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

Monitors



Oximeters

GE Ohmeda TruSat

Accessories

OXY-BC-5	Bedside clip (5/box)
OXY-HB	Replacement headband
OXY-RTW	Wide adhesive replacement tape for use with Wrap, Sensitive Skin, or AllFit sensors
OXY-RWL	Foam wrap replacement, large, weight range ≥ 3 kg (6.6 lb) for use with Wrap sensor
OXY-RWM	Foam wrap replacement, large, weight range ≥ 3 kg (6.6 lb) for use with Sensitive Skin sensor
OXY-RWS	Foam wrap replacement, small, weight range < 3 kg (6.6 lb) or small finger, for use with Sensitive Skin sensor
OXY-RTB	Replacement adhesive tape (bears), for use with AllFit sensor
OXY-RT	Replacement adhesive tape (blue), for use with AllFit sensor
OXY-SND	Infant foam sandal, for use with Sensitive Skin sensor

Connecting cables

TruSat connector



Interconnect cable with TruSat connector (1/box)

TS-M3 3 m (10 ft)

Sensor type

Reusable sensors

TS-E-D	TruSignal Ear sensor 1 m (3.3 ft) (1/box), requires interconnect cable
TS-F-D	TruSignal Finger sensor 1 m (3.3 ft) (1/box), requires interconnect cable
TS-W-D	TruSignal Wrap sensor 1 m (3.3 ft) (1/box), requires interconnect cable
TS-SE-3	TruSignal Sensitive Skin sensor 1 m (3.3 ft) (3/box), requires interconnect cable
TS-SA-D	TruSignal FingerTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable
TS-SP-D	TruSignal PediTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable

Adhesive sensors

TS-AF-10 TS-AF-25	TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable TS-AF-10 - 10/box, TS-AF-25 - 25/box
TS-AP-10 TS-AP-25	TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable TS-AP-10 - 10/box, TS-AP-25 - 25/box
TS-AAW-10 TS-AAW-25	TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable TS-AAW-25 - 10/box, TS-AAW-25 - 25/box
TS-PAW-10 TS-PAW-25	TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable TS-PAW-10 - 10/box, TS-PAW-25 - 25/box

Monitors



1. Oximeters

GE Ohmeda TuffSat
3775, 3700 Series, 3800, 3900 and 3900P

2. Configured

Cardiacap/5 with N-XOSAT option

3. Modular monitors with modules

S/5, AS/3, CS/3 modular monitors with
M-OSAT module

Accessories

OXY-BC-5	Bedside clip (5/box)
OXY-HB	Replacement headband
OXY-RTW	Wide adhesive replacement tape for use with Wrap, Sensitive Skin, or AllFit sensors
OXY-RWL	Foam wrap replacement, large, weight range ≥ 3 kg (6.6 lb) for use with Wrap sensor
OXY-RWM	Foam wrap replacement, large, weight range ≥ 3 kg (6.6 lb) for use with Sensitive Skin sensor
OXY-RWS	Foam wrap replacement, small, weight range < 3 kg (6.6 lb) or small finger, for use with Sensitive Skin sensor
OXY-RTB	Replacement adhesive tape (bears), for use with AllFit sensor
OXY-RT	Replacement adhesive tape (blue), for use with AllFit sensor
OXY-SND	Infant foam sandal, for use with Sensitive Skin sensor

Connecting cables

Ohmeda connector



Interconnect cable with Ohmeda connector (1/box)

TS-H3 3 m (10 ft)



No interconnect cable required

Sensor type

Reusable sensors

- TS-F-D** TruSignal Finger sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-W-D** TruSignal Wrap sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-SE-3** TruSignal Sensitive Skin sensor 1 m (3.3 ft) (3/box), requires interconnect cable
- TS-SA-D** TruSignal FingerTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-SP-D** TruSignal PediTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable

Adhesive sensors

- TS-AF-10** TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable
- TS-AF-25** TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable
- TS-AP-10** TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable
- TS-AP-25** TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable
- TS-AAW-10** TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable
- TS-AAW-25** TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable
- TS-PAW-10** TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable
- TS-PAW-25** TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable

Integrated sensors

- TS-F1-H** Integrated Finger sensor with Ohmeda connector, 1 m (3.3 ft), (1/box)
- TS-E4-H** Integrated Ear sensor with Ohmeda connector, 4 m (13 ft), (1/box) headband and clip included

Monitors



1. Configured

Cardiacap/5
S/5 Light Monitor

2. Modular monitors with modules

CARESCAPE Monitor B850 with M-PRESTN..01 or M-RESTN..01 modules, S/5, AS/3, CS/3 modular monitors with hemodynamic M-modules

Accessories

OXY-BC-5	Bedside clip (5/box)
OXY-HB	Replacement headband
OXY-RTW	Wide adhesive replacement tape for use with Wrap, Sensitive Skin, or AllFit sensors
OXY-RWL	Foam wrap replacement, large, weight range ≥ 3 kg (6.6 lb) for use with Wrap sensor
OXY-RWM	Foam wrap replacement, large, weight range ≥ 3 kg (6.6 lb) for use with Sensitive Skin sensor
OXY-RWS	Foam wrap replacement, small, weight range < 3 kg (6.6 lb) or small finger, for use with Sensitive Skin sensor
OXY-RTB	Replacement adhesive tape (bears), for use with AllFit sensor
OXY-RT	Replacement adhesive tape (blue), for use with AllFit sensor
OXY-SND	Infant foam sandal, for use with Sensitive Skin sensor

Connecting cables

Datex connector



Interconnect cable with Datex connector (1/box)

TS-N3 3 m (10 ft)



No interconnect cable required

Sensor type

Reusable sensors

- TS-F-D** TruSignal Finger sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-W-D** TruSignal Wrap sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-SE-3** TruSignal Sensitive Skin sensor 1 m (3.3 ft) (3/box), requires interconnect cable
- TS-SA-D** TruSignal FingerTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-SP-D** TruSignal PediTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable

Adhesive sensors

- TS-AF-10** TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable
- TS-AF-25** TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable
- TS-AP-10** TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable
- TS-AP-25** TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable
- TS-AAW-10** TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable
- TS-AAW-25** TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable
- TS-PAW-10** TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable
- TS-PAW-25** TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable

Integrated sensors

- TS-F4-H** Integrated Finger sensor with Datex connector, 1 m (3.3 ft), (1/box)
- TS-E4-N** Integrated Ear sensor with Datex connector, 4 m (13 ft), (1/box) headband and clip included

Monitors



1. Configured

CARESCAPE V100 with GE TruSignal option
 CARESCAPE VC150 with GE TruSignal option
 B30 monitor

2. Modular monitors with modules

CARESCAPE modular monitors with E-PRESTN,
 E-RESTN, E-PSM(P) module, CARESCAPE ONE with
 TruSignal Option

FM and FM Light with E-PSM(P) module

S/5 modular monitors with E-PRESTN, E-RESTN,
 E-PSM(P) module

3. PROCARE monitors

PROCARE* Monitor B20(i)

PROCARE Monitor B40(i)

Accessories

- OXY-BC-5** Bedside clip (5/box)
- OXY-HB** Replacement headband
- OXY-RTW** Wide adhesive replacement tape for use with Wrap, Sensitive Skin, or AllFit sensors
- OXY-RWL** Foam wrap replacement, large, weight range ≥ 3 kg/6.6 lb for use with Wrap sensor
- OXY-RWM** Foam wrap replacement, large, weight range ≥ 3 kg/6.6 lb for use with Sensitive Skin sensor
- OXY-RWS** Foam wrap replacement, small, weight range < 3 kg/6.6 lb or small finger, for use with Sensitive Skin sensor
- OXY-RTB** Replacement adhesive tape (bears), for use with AllFit sensor
- OXY-RT** Replacement adhesive tape (blue), for use with AllFit sensor
- OXY-SND** Infant foam sandal, for use with Sensitive Skin sensor

Connecting cables

GE connector



Interconnect cable with GE connector (1/box)

TS-G3 3 m (10 ft)

CARESCAPE ONE connector



Interconnect cable with CARESCAPE ONE connector (1/box)

2097648-001 1.8 m (6 ft)

2097648-002 3 m (10 ft)

GE connector



No interconnect cable required

Sensor type

Reusable sensors

- TS-E-D** TruSignal Ear sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-F-D** TruSignal Finger sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-W-D** TruSignal Wrap sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-SE-3** TruSignal Sensitive Skin sensor 1 m (3.3 ft) (3/box), requires interconnect cable
- TS-SA-D** TruSignal FingerTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable
- TS-SP-D** TruSignal PediTip sensor 1 m (3.3 ft) (1/box), requires interconnect cable

Adhesive sensors

- TS-AF-10** TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable
- TS-AF-25** TruSignal AllFit sensor 0.5 m (1.6 ft), extra adhesive tape included, requires interconnect cable
- TS-AP-10** TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable
- TS-AP-25** TruSignal Adult/Pediatric sensor 0.3 m (1ft), requires interconnect cable
- TS-AAW-10** TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable
- TS-AAW-25** TruSignal Adult sensor 0.3 m (1 ft), requires interconnect cable
- TS-PAW-10** TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable
- TS-PAW-25** TruSignal Pediatric/Infant sensor 0.3 m (1 ft), requires interconnect cable

Integrated sensors

- TS-E2-GE** TruSignal Integrated Ear sensor with GE connector 2 m (6.6 ft), (1/box), headband and clip included
- TS-E4-GE** TruSignal Integrated Ear sensor with GE connector 4 m (13 ft), (1/box), headband and clip included
- TS-F2-GE** TruSignal Integrated Finger sensor with GE connector 2 m (6.6 ft), (1/box)
- TS-F4-GE** TruSignal Integrated Finger sensor with GE connector 4 m (13 ft), (1/box)
- TS-SA4-GE** TruSignal Integrated FingerTip sensor with GE connector 4 m (10 ft), (1/box)
- TS-SP3-GE** TruSignal Integrated PediTip sensor with GE connector 3 m (10 ft), (1/box)


Monitors

1. Configured

CARESCAPE V100 with Masimo option
 CARESCAPE VC150 with Masimo option
 Corometrics 250 Series with Masimo option
 B20(i) with Masimo option
 B40(i) with Masimo option
 B105, B125 with Masimo option
 Dash 5000 with Masimo option

2. Modular monitors with modules

CARESCAPE Monitor B850 with CARESCAPE Patient Data Module with Masimo option with E-MASIMO module
 CARESCAPE Monitor B650 with CARESCAPE Patient Data Module with Masimo option with E-MASIMO module
 CARESCAPE Monitor B450 with CARESCAPE Patient Data Module with Masimo option with E-MASIMO module
 Solar 8000i v5 with CARESCAPE Patient Data Module with Masimo option



Accessories

MASIMO-RD 4201

RD SET™ Training Kit. Includes: 1 RD SET Adt, 1 RD SET Pdt, 1 RD SET Inf, 1 RD SET Neo, 1 RD SET NeoPt, 1 neonatal foot with application card, 1 RD SET sensor application card, Direction for use

MASIMO-RD 4046

RD SET Inf Replacement Tapes. For use with RD SET Inf. 102/box, non-sterile. Not made with natural rubber latex.

MASIMO-RD 4047

RD SET Neo Replacement Tapes. For use with RD SET Neo. 102/box, non-sterile. Not made with natural rubber latex.

MASIMO-RD 4048

RD SET NeoPt Replacement Wraps. For use with RD SET NeoPt. 10/pack, non-sterile. Not made with natural rubber latex.

Connecting cables


GE connector



Patient cable with GE connector

MASIMO-RD 4084 1.5 m (5 ft)
MASIMO-RD 4085 3.7 m (12 ft)

VC150 connector



Patient Cable with RD Rainbow MD20 connector

MASIMO-RD 4103 1.5 m (5 ft)
MASIMO-RD 4104 3.7 m (12 ft)

Adapter cables



Patient Cable with LNOP® PC connector

MASIMO-RD 4088 0.5 m (1.5 ft)



Patient Cable with LNC connector


MASIMO-RD 4105 0.5 m (1.5 ft)



Patient Cable with M-LNC connector

MASIMO-RD 4090 0.5 m (1.5 ft)

Sensor type



Reusable sensors

MASIMO-RD 4050 RD SET DCI Adult Reusable Sensor, 3 ft. 1/box, non-sterile

MASIMO-RD 4051 RD SET DCI-P Pediatric/Slender Digit Reusable Sensor, 3 ft. 1/box, non-sterile.

MASIMO-RD 4053 RD SET TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, non-sterile.

Adhesive sensors

MASIMO-RD 4000 RD SET Adt Adult SpO₂ Adhesive Sensor. Single-patient-use. 20/box, non-sterile.

MASIMO-RD 4001 RD SET Pdt Pediatric SpO₂ Adhesive Sensor. Single-patient-use. 20/box, non-sterile

MASIMO-RD 4002 RD SET Inf Infant SpO₂ Adhesive Sensor. Single-patient-use. 20/box, non-sterile.

MASIMO-RD 4003 RD SET Neo Adult/Neonatal SpO₂ Adhesive Sensor. Single-patient-use. 20/box, non-sterile.

MASIMO-RD 4004 RD SET NeoPt Neonatal SpO₂ Adhesive Sensor. Single-patient-use. 20/box, non-sterile.

MASIMO-RD 4005 RD SET NeoPt-500 Neonatal SpO₂ Non-adhesive Sensor. Single-patient-use. 20/box, non-sterile.

Monitors



1. Configured

CARESCAPE V100 with Masimo option
 CARESCAPE VC150 with Masimo option
 Corometrics 250 Series with Masimo option
 Dash 2500 with Masimo option
 Dash 3000/4000/5000 with Masimo option

2. Modular monitors with modules

S/5, AS/3, CS/3 modular monitors with E-MASIMO module
 CardiLab/Mac-Lab with Tram module 451M
 CARESCAPE Monitor B850
 with CARESCAPE ONE with Masimo option
 with CARESCAPE Patient Data Module with Masimo option
 with E-MASIMO module
 with Masimo SpO2 module for Tram-Rac
 with Tram module 451M/851M
 CARESCAPE Monitor B650
 with CARESCAPE ONE with Masimo option
 with CARESCAPE Patient Data Module with Masimo option
 with E-MASIMO module
 CARESCAPE Monitor B450
 with CARESCAPE ONE with Masimo option
 with CARESCAPE Patient Data Module with Masimo option
 Solar 8000M v5/8000i v5
 with CARESCAPE Patient Data Module with Masimo option
 with Tram module 451M/851M
 Solar 8000/8000M/8000i
 with Masimo SpO2 module for Tram-Rac
 with Tram module 451M/851M
 Solar 9500 with Tram module 451M/851M

3. PROCARE monitors

PROCARE Monitor B20(i) with Masimo option
 PROCARE Monitor B40(i) with Masimo option

Accessories

- 2051071-001** Replacement tape for LNCS Inf series sensors (100/box)
- 2051073-001** Replacement tape for Neo series sensors (100/box)
- 2051074-001** Replacement tape for NeoPt series sensors (10/box)

Connecting cables

GE connector



Masimo LNC-10 interconnect cable for Dash, Tram, Corometrics

2027263-002 3 m (10 ft)

CARESCAPE ONE connector




Interconnect cable with CARESCAPE ONE connector (1/box)

2103987-001 1.8 m (6 ft)

2103987-002 0.9 m (3 ft)

Sensor type



Reusable sensors

2027258-001 Masimo LNCS DCI Reusable sensor, Adult

2027259-001 Masimo LNCS DCIP Reusable sensor, Pediatric

2027261-001 Masimo LNCS TC-I TipClip Reusable Ear sensor

Adhesive sensors

2027253-001 Masimo LNCS Adhesive sensor, LNCS-Adtx, Adult (20/box)

2027254-001 Masimo LNCS Adhesive sensor, LNCS-Pdtx, Pediatric (20/box)

2051068-001 Masimo LNCS Adhesive sensor, LNCS-Inf, Infant (20/box)

2051069-001 Masimo LNCS Adhesive sensor, LNCS-Neo Neonatal (20/box)

2051070-001 Masimo LNCS Adhesive sensor, LNCS-Neo-Pt, Neonatal (20/box)

Monitors



Configured

- GE DINAMAP ProCare with Masimo option
- GE DINAMAP PRO Series v2 with Masimo option
- GE DINAMAP PRO 1000 v3 with Masimo option

Accessories

- 2051071-001** Replacement tape for LNCS Inf series sensors (100/box)
- 2051073-001** Replacement tape for Neo series sensors (100/box)
- 2051074-001** Replacement tape for NeoPt series sensors (10/box)


Connecting cables



Masimo LNC-10 interconnect cable for GE DINAMAP

2027263-001 3 m (10 ft)

Sensor type



Reusable sensors

- 2027258-001** Masimo LNCS DCI Reusable sensor, Adult
- 2027259-001** Masimo LNCS DCIP Reusable sensor, Pediatric
- 2027261-001** Masimo LNCS TC-I TipClip Reusable Ear sensor

Adhesive sensors

- 2027253-001** Masimo LNCS Adhesive sensor, LNCS-Adtx, Adult (20/box)
- 2027254-001** Masimo LNCS Adhesive sensor, LNCS-Pdtx, Pediatric (20/box)
- 2051068-001** Masimo LNCS Adhesive sensor, LNCS Inf, Infant (20/box)
- 2051069-001** Masimo LNCS Adhesive sensor, LNCS Neo Neonatal (20/box)
- 2051070-001** Masimo LNCS Adhesive sensor, LNCS-Neo-Pt, Neonatal (20/box)

Monitors



OxiMax

CARESCAPE Monitor B850
with CARESCAPE Patient Data Module with Nellcor OxiMax™ option
with Tram module 451N5/851N5
with CARESCAPE ONE with Nellcor Option

CARESCAPE Monitor B650
with CARESCAPE Patient Data Module with Nellcor OxiMax option
with CARESCAPE ONE with Nellcor Option

CARESCAPE Monitor B450
with CARESCAPE Patient Data Module with Nellcor OxiMax option
with CARESCAPE ONE with Nellcor Option

CARESCAPE V100 with Nellcor OxiMax option
CARESCAPE VC150 with Nellcor OxiMax option

Cardiolab/Mac-Lab with Tram module 451N5

Dash 2500 with Nellcor OxiMax option
Dash 3000/4000/5000 with Nellcor OxiMax option

Solar 8000/8000M/8000i and 9500 with Tram module 451N5/851N5

Solar 8000M v5/8000i v5 with CARESCAPE Patient Data Module with Nellcor OxiMax option

E-NSATX module for S/5, AS/3, CS/3 modular monitors and CARESCAPE modular monitors

OxiMax



GE DINAMAP ProCare with Nellcor OxiMax option

GE DINAMAP PRO Series v2 with Nellcor OxiMax option

GE DINAMAP PRO 1000 v3 with Nellcor OxiMax option

E-NSAT module for S/5, AS/3, CS/3 modular monitors and CARESCAPE modular monitors

Corometrics 250 Series
with Nellcor OxiMax option



Connecting cables

GE connector



Interconnect cable with GE connector and Nellcor OxiMax connector

2021406-001 3 m (10 ft)

2021406-002 1.2 m (4 ft)

CARESCAPE ONE connector



Interconnect cable with CARESCAPE ONE connector (1/box)

2082747-001 1.2 m (4 ft)

2082747-002 3.6 m (12 ft)



Interface cable Nellcor OxiMax DOC-10

2008773-001 2.9 m (9.5 ft)


GE connector



Interconnect cable with GE connector and Nellcor OxiMax connector (for Corometrics only)

2025350-001 3 m (10 ft)

Sensor type



Reusable sensors

70124021/407705-006 Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)

70124033 Nellcor SpO₂ Multisite Sensor D-YS

70124034 Nellcor SpO₂ Ear-Clip D-YSE, for use with 70124033 (1/box)

70124035 Nellcor SpO₂ Sensor OXY-A/N, Neonatal/Adult (1/box)

414248-001 Nellcor SpO₂ Sensor OXI-P/I, Pediatric/Infant (1/box)

Adhesive sensors

70124022 Nellcor SpO₂ Sensor MAX-P, Adhesive, Pediatric (24/box)

70124026 Nellcor SpO₂ Sensor MAX-I, Adhesive, Infant (24/box)

70124027 Nellcor SpO₂ Sensor MAX-A, Adhesive, Adult (24/box)

70124032 Nellcor SpO₂ Sensor MAX-N Adhesive, Neonatal (24/box)

407705-005 Nellcor SpO₂ Sensor MAX-R, Adhesive, Nasal (24/box)

2028117-001 Nellcor SpO₂ sensor MAX-AL, Adhesive, Adult Long (24/box)

Monitors



OxiSmart

CARESCAPE Monitor B850 with Tram module¹
451/451N/851/851N

CardioLab/Mac-Lab with Tram module 451/451N

Corometrics 120/2120 Series

Dash 3000/4000 (v5 or lower) with legacy Marquette option

Dash 3000/4000/5000 with GE Ohmeda option

Responder 3000 with Nellcor OxiSmart option

Solar 8000/8000M/8000i and 9500 with Tram module 451/451N/851/851N



M-NSAT module for S/5, AS/3, CS/3 modular monitors and CARESCAPE modular monitors[†]

Configured monitor:
Cardiacap/5 with N-XNSAT Nellcor option

Corometrics 118

[†] Refer to the CARESCAPE Monitor B850 and CARESCAPE Monitor B650 User Manuals for compatible M-NSAT module versions.


Connecting cables

GE connector



Interconnect cable with GE connector and Nellcor OxiSmart connector

2006644-001 2.9 m (9.5 ft)



Interconnect cable Nellcor MC-10


903131 3 m (10 ft)



Interconnect cable Nellcor Corometrics 118 monitor

4033CA 3 m (10 ft)

Sensor type



Reusable sensors

70124021/407705-006 Nellcor DuraSensor Reusable Finger Probe DS100A (1/box)

70124033 Nellcor SpO₂ Multisite Sensor D-YS

70124034 Nellcor SpO₂ Ear-Clip D-YSE, for use with 70124033 (1/box)

70124035 Nellcor SpO₂ Sensor OXY-A/N, Neonatal/Adult (1/box)

414248-001 Nellcor SpO₂ Sensor OXI-P/I, Pediatric/Infant (1/box)

Adhesive sensors

70124022 Nellcor SpO₂ Sensor MAX-P, Adhesive, Pediatric (24/box)

70124026 Nellcor SpO₂ Sensor MAX-I, Adhesive, Infant (24/box)

70124027 Nellcor SpO₂ Sensor MAX-A, Adhesive, Adult (24/box)

70124032 Nellcor SpO₂ Sensor MAX-N Adhesive, Neonatal (24/box)

407705-005 Nellcor SpO₂ Sensor MAX-R, Adhesive, Nasal (24/box)

2028117-001 Nellcor SpO₂ sensor MAX-AL, Adhesive, Adult Long (24/box)

Note: This guide contains SpO₂ cables and sensors for different technologies. Always refer to the appropriate monitoring device's User Guide for the accurate list of compatible cables and sensors

✦ All Nellcor entries in this guide are subject to limited availability. Please consult your GE representative.

¹ Not compatible with software version 3

Monitors



CARESCAPE Telemetry T4 Transmitter
ApexPro CH Telemetry Transmitter
CARESCAPE Telemetry T14 Transmitter



ApexPro FH Telemetry Transceiver

Connecting cables



Interconnect cable for CH, T4 & T14:

2062190-001 uSpO₂ Pulse Oximetry Cable for GE ApexPro, GE ApexPro CH, CARESCAPE Telemetry



Interconnect cable for FH:

2062190-002 uSpO₂ Pulse Oximetry Cable For GE ApexPro FH

Sensor type

Reusable sensors

- 2066338-001** LNCS DB-I Adult Soft Reusable Sensor (30 kg and greater)
- 2066339-001** LNCS Y-I Multisite Reusable Sensor
- 2027258-001** Sensor LNCS DC-1 Adult Reusable
- 2027259-001** Sensor LNCS DC-1P PED Reusable
- 2027261-001** Sensor LNCS TC-1 Tip-Clip Ear Reusable

Adhesive sensors

- 2066333-001** LNCS[®] Adtx-3 Adult SpO₂ Adhesive Sensors (30 kg and greater)
- 2066334-001** LNCS Pdtx-3 Pediatric SpO₂ Adhesive Sensors (10 kg - 50 kg)
- 2066335-001** LNCS Inf-3 Infant SpO₂ Adhesive Sensors (3 kg - 20 kg)
- 2066336-001** LNCS Neo-3 SpO₂ Adhesive Sensors (40 kg and greater)

Note: This guide contains SpO₂ cables and sensors for different technologies. Always refer to the appropriate monitoring device's User Guide for the accurate list of compatible cables and sensors

♣ All Nellcor entries in this guide are subject to limited availability. Please consult your GE representative.

† Refer to the CARESCAPE Monitor B850 and CARESCAPE Monitor B650 User Manuals for compatible M-NSAT module versions.

GE Healthcare is a leading global medical technology and digital solutions innovator. GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world. Follow us on Facebook, LinkedIn, Twitter and Insights, or visit our website www.gehealthcare.co.uk for more information.

© 2021 General Electric Company - All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE, GE Monogram, ApexPro, Cardiocap, CARESCAPE, Corometrics, DINAMAP, Marquette, Ohmeda, Solar, TruSignal, TruSat and TuffSat are trademarks of General Electric Company.

MASIMO is a trademark of Masimo Corporation.

Nellcor, OxiMax and OxiSmart are trademarks of a Medtronic Company.

JB01066XE 01/2021