


TOMORROW TODAY



AIR™

Simply better

gehealthcare.com/mr





Simply better

It's time for a simply better MR experience. An experience that embraces the clinical power of MR with a more ubiquitous design. For us, this is AIR™.

AIR™ is more than a name. It's a way of thinking. A symbol for how we approach both the form and the function of MR. From coils to workflow to image quality, we deconstruct everything down to its most fundamental, simplified version. And improve it. Each new innovation builds off of this core philosophy, increasing the value of the greater AIR™ ecosystem.

Holistically, AIR™ is MR. Simply better.



Welcome to the AIR™ family


The AIR™ family of products delivers clinical versatility and comfort, intelligent productivity improvements and consistently superior image quality.



AIR™ Coils

It started with an engineering breakthrough that led to the industry's first truly lightweight coil design. A design that conforms to the human body like a comforting blanket completely changing the patient experience. But this was just the start.

■ Clinical versatility and comfort



AIR™ Workflow

The technology at the core of this new coil design laid the groundwork for intelligent applications that automate and personalize fundamental steps in the MR workflow, so you can expect consistent image quality in predictable scan times from any technologist.


■ Intelligent productivity improvements



AIR™ Image Quality

Now, we're adding reconstruction software that reduces background noise and out-of-FOV artifacts to foundationally improve the quality of your MR images.

■ Consistently superior image quality



AIR™ Coils

Clinical versatility and comfort


Awarded Best New Radiology Device of 2019, AIR™ Coils are the foundation of a simply better MR experience. The engineering breakthrough at the heart of our AIR™ Coils allowed us to create a revolutionary coil design that is lighter, offers more flexibility and provides greater coverage, laying the groundwork for greater positioning freedom and a comfortable patient experience.



Lightweight and bendable, this series of linked resonators replaces what used to be rigid circuit boards and lumped components. It's the flexible core of our form-fitting coil design.



Neck Imaging
AIR™ 48ch Head Coil + AIR™ Anterior Array Coil +
Posterior Array, 3D T1 HyperCube FLEX with HS
1 x 0.8 x 1.2 mm, 4:21 min



AIR™ 48ch Head Coil*

Designed in collaboration with the NFL** for high-performance brain imaging.

- Adjustable to fit 99.99 percent of patients
- Fit-adaptable, comfort tilt design
- Improved parallel imaging





* For 3.0T DLI systems.
** In 2013, GE and the National Football League (NFL) teamed up to launch the Head Health Initiative, a collaboration to accelerate diagnosis and improve treatment for traumatic brain injury.



Automatic coil element selection for the prostate with AIR Touch™

Axial T2 FSE
0.7 x 0.7 x 2.5 mm




AIR™ Anterior Array Coil

Scan the chest, abdomen and pelvis without repositioning the coil.

- 65 cm of coverage
- 30-channel design
- Weighs less than 0.35 grams per cm²






AIR™ Workflow

Intelligent productivity improvements

Enhance your MR productivity with intelligent workflow applications developed to optimize your scans with accelerated scan times, increase diagnostic confidence across skill levels and consistently deliver accurate results. Automated applications, AIR Touch™ and AIR x™, make a clinically impactful difference for a simply better workflow.






AIR Touch™

Smart coil element selection automatically knows the best combination for every patient.

- Simplifies scan setup with customized parameters
- Optimizes uniformity and SNR
- Automatically selects which coils and elements to use



AIR Touch™ acts as the bridge between the coils and the system, allowing you to landmark your patient with a single touch.




AIR x™

Intelligent MR slice prescription for routine and challenging neurological exams.

- Powered by a deep learning algorithm created from a database of 36,000 images
- Automatically detects anatomy and prescribes slices in the brain
- Delivers consistent and quantifiable results
- Helps eliminate rescans and scanning inefficiencies




Powered By
Edison



Time point one

Time point two

Time point three



AIR™ Image Quality

Consistently superior image quality

AIR™ Image Quality completes the AIR™ family of products with image reconstruction software that reduces background noise and out-of-FOV artifacts. It helps improve signal-to-noise in every image without having to overcompensate in your scanning protocol.

AIR™ Recon

Makes exceptional image quality in faster scan times the new standard for MR imaging.

- Available on the majority of routine clinical scans
- Seamless to use and always on
- Reduces background noise and out-of-FOV artifacts
- Improves signal-to-noise ratio without having to overcompensate in your protocols



Conventional



AIR™ Recon

Clearly consistent images with AIR™ Recon




Conventional



AIR™ Recon

AIR™ 48ch Head Coil
0.5 x 0.5 x 2.5 mm
3:56 min



Conventional



AIR™ Recon

AIR™ Anterior Array Coil
0.5 x 0.5 x 3 mm
5:02 min

Get closer to what you need to see

AIR™ combined with the advanced applications of SIGNA™Works enables diagnostic confidence in every scan through industry-leading image quality and imaging capabilities.

AIR™ Anterior Array Coil



Two-station chest/abdomen/pelvis acquired with only one AIR Anterior Array coil without moving the coil or subject




Sagittal T2 PROPELLER MB
0.7 x 0.7 x 3.5 mm




Axial Inhance IFIR
Non-contrast
Coronal reformat
0.8 x 1.3 x 1.6 mm


AIR™ 48ch Head Coil




Volume rendering
HyperBand DTI FiberTrak fused with HyperSense 3D TOF



HyperSense
3D TOF
0.3 x 0.3 x 0.5 mm



3D BRAVO
1 x 1 x 1 mm



MUSE DWI
b1000
0.9 x 0.8 x 5 mm

Coil combinations with AIR Touch™


AIR Touch™ automatically selects the best coil combination and utilizes anatomical-based protocol optimization to personalize care and improve exam times by 53 percent.*



AIR™ 48ch Head Coil +
AIR™ Anterior Array Coil +
Posterior Array
Sagittal T2 frFSE
0.6 x 0.9 x 3 mm



AIR™ Anterior Array +
Posterior Array
Axial T2 PROPELLER
0.6 x 0.6 x 3 mm



AIR™ 48ch Head Coil +
AIR™ Anterior Array Coil +
Posterior Array
Coronal 3D T2 STIR
HyperCube with HyperSense
1.25 x 1.25 x 1.6 mm





Curved reformat Sagittal reformat
Axial reformat

Two AIR™ Anterior Array Coils +
Posterior Array
Coronal STIR FSE
1 x 1.8 x 2.5 mm

*Compared to previous workflows without AIR Touch™.

It all links up to a simply better MR experience

What really sets our AIR™ family of products apart is the way they all work together. It started with an engineering breakthrough that led to the industry's first truly lightweight coil design and laid the groundwork for intelligent workflow applications and an all-new reconstruction algorithm. Each one builds off of the potential of the rest of the family, creating a chain of innovation that all links together to transform the entire MR experience.



AIR™
Simply better



For more information, visit gehealthcare.com/AIR
or contact your GE Healthcare Sales Representative.

GE Healthcare is a leading provider of medical imaging, monitoring, biomanufacturing, and cell and gene therapy technologies. GE Healthcare enables precision health in diagnostics, therapeutics and monitoring through intelligent devices, data analytics, applications and services. With over 100 years of experience and leadership in the healthcare industry and more than 50,000 employees globally, GE Healthcare helps healthcare providers, researchers and life sciences companies in their mission to improve outcomes for patients around the world. Follow us on Facebook, LinkedIn, Twitter and The Pulse for latest news, or visit our website www.gehealthcare.com for more information.

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, AIR, SIGNA, AIR Touch, AIR x, Imagination at work and the GE Monogram are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB76849FR