



Clinical accessories portfolio

MAC VU360

ECG Accessories with IEC standard

ECG CABLES AND LEADWIRES



PRODUCT	SIZE	PART NUMBER	QUANTITY
ECG Trunk Cable CC14	2,5 m	2096724-067	1
Leadwire Set, Base 10, banana, IEC Limb leads 125 cm/ Chest leads 80 cm	mixed length	2104724-002	1
Leadwire Set, 4 Add-on, banana, IEC-HEIM	80 cm	2104750-001	1
Leadwire Set, Base 10, grabber, IEC Limb leads 125 cm/ Chest leads 80 cm	mixed length	2104748-001	1
Leadwire Set, 4 Add-on, grabber, IEC-HEIM	80 cm	2104725-002	1



ACCESSORIES

PRODUCT	PART NUMBER	QUANTITY
Adapter Clip, Wide Mouth Clip	2104729-001	14
ECG Leadwire port plug	2104767-001	4

The list of accessories displayed in this catalogue is not exhaustive, Contact your GE Healthcare representative to obtain the list of all accessories compatible with your equipment.



PAPER

PRODUCT	PART NUMBER	QUANTITY
Thermal paper A4, red grid 205mm wide, z-fold, block queue, 150 sheets	2104772-001	10



ELECTRODES

PRODUCT	FORM	SIZE	PART NUMBER	QUANTITY
ECG Snap Adult Electrode (tape) - Solid Hydrogel	Teardrop	47,625 mm	2104777-001	600
ECG Snap Adult Electrode (foam) - 10% KCl Wet gel	Teardrop	38.1 mm	2104779-001	600
ECG Tab Adult Electrode (Adhesive) - Solid Hydrogel	Medium Tab	34.036 mm x 22.86 mm	2106911-001	4000
ECG Tab Adult Electrode (Adhesive) - Solid Hydrogel	Large Tab	28.575 mm x 36.525 mm	2104780-001	4000
ECG Tab Pediatric Electrode (Adhesive) - Adhesive Hydrogel	Small Tab	12.7 mm x 32.54 mm	2104778-001	4000

GE imagination at work

About GE Healthcare

GE Healthcare is a leading provider of medical imaging, monitoring, biomanufacturing, and cell and gene therapy technologies. GE Healthcare enables precision health in diagnostics, therapeutics and monitoring through intelligent devices, data analytics, applications and services. With over 100 years of experience and leadership in the healthcare industry and more than 50,000 employees globally, GE Healthcare helps healthcare providers, researchers and life sciences companies in their mission to improve outcomes for patients around the world. Follow us on Facebook, LinkedIn, Twitter and The Pulse for latest news, or visit our website www.gehealthcare.com for more information.

The list of accessories displayed in this catalogue is not exhaustive, Contact your GE Healthcare representative to obtain the list of all accessories compatible with your equipment.

www.gehealthcare.com.

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information.

Data subject to change.

GE, the GE Monogram and Imagination at work are trademarks of General Electric Company.

*All other third-party trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE.

Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

© 2019 General Electric Company. All rights reserved.

JB65724XE