



# Easy-Fil Vaporizer Bottle Adapters

Easy-Fil vaporizer bottle adapters from GE Healthcare are tested and verified to function with compatible GE Healthcare vaporizers. The adapters perform optimally in filling vaporizers with the associated anesthetic agent while also offering the highest strength, quality and durability product.

- **Strength:** New stringent testing completed to ensure part consistency for the highest strength possible
- **Quality:** All parts are 100% leak tested with three different methods to ensure bottle to vaporizer seal
- **Durability:** Tested to over 1,000 fill cycles while still compliant with all other testing to ensure long lifespan

## Product compatibility

- Easy-Fil Bottle Adapter, Vaporizer, Isoflurane: All isoflurane vaporizers with Easy-Fil filler
- Easy-Fil Bottle Adapter, Vaporizer, Halothane: All halothane vaporizers with Easy-Fil filler
- Easy-Fil Bottle Adapter, Vaporizer, Sevoflurane: All sevoflurane vaporizers with Easy-Fil filler

Our equipment.  
Our accessories.  
Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

## PRODUCT SPECIFICATIONS

### Material composition

---

(All Easy-Fil bottle adapters are latex-free)

Body and Retainer	Nylon 6/6 with glass fill with colorant
Stem	LDPE
Spring	Stainless Steel
O-Ring	Teflon

### Suggested lifespan

---

1,000 fill cycles

### Storage Requirements

---

No specific or required storage criteria

### ISO Standards

---

Parts tested to ISO 5360:2016. Additional Internal testing performed to surpass ISO requirements.

## PART LIST

Part Number	Description
2106572-001	Easy-Fil Bottle Adapter, Vaporizer, Isoflurane
2106572-002	Easy-Fil Bottle Adapter, Vaporizer, Halothane
2106572-003	Easy-Fil Bottle Adapter, Vaporizer, Sevoflurane

**Product may not be available in all countries and regions.** Contact a GE Healthcare Representative for more information. Please visit [www.gehealthcare.com](http://www.gehealthcare.com).

*Distributed by:*

Datex-Ohmeda, Inc.



*Manufactured by:*

Datex-Ohmeda, Inc.  
3030 Ohmeda Drive  
Madison, WI 53718

## Imagination at work

© 2019 General Electric Company - All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Imagination at work are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2192541 1/19